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Reusable tourniquet

Device description

The device immobilises the blood flow in peripheral venous vessels to expand the immobilised blood vessel and enable improved vein puncture. The device can also fix the vessel for puncture without immobilising the blood flow.

This product is a medical device. The device may only be used by trained medical personnel.

Patient group

The device is intended for use on adults and children.

Materials

Aluminium, cotton, PA, PUR

Contraindications

The use of the device is contraindicated in patients with known sensitivities or allergies to the above materials.

The device should be used only on intact skin that is not affected by disease.

Application

Fig. 1: Place the tourniquet around the limb several centimetres proximal from the puncture site. Insert the closure into the lock until it clicks shut.

Fig. 2: Tighten the tourniquet and push the lock to the extremity to immobilise the blood flow in the peripheral vein.

Fig. 3: The opening in the closure can be placed over the vein that is to be punctured to fix it in position for puncturing.

The blood vessel in this case is not immobilised and an injection can be administered with the tourniquet in place.

Fig. 4: To loosen the tourniquet or to release the immobilisation gradually press the metal button on the lock.

Fig. 5: For complete release press the closure to unlock it and release it from the lock.

Directions

● After use, the device must be disinfected with standard chemicals.

● Dispose of the device in accordance with your internal instructions and local regulations.

● Do not use the tourniquet for bandaging of bleeding, e.g. from traumatic injuries or failed punctures.

● When placing the tourniquet be careful not to trap or injure skin or skin folds.

● Do not allow the closure to unlock suddenly, because this may cause injury.

Additional information

All serious incidents occurring in connection with the medical device must be reported to the manufacturer and the responsible authorities.

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