# FLOQSwabs®

Package Insert and How to use guide  
for the medical laboratory

CE 0123

**FLOQSwabs® for self-collection of vaginal swab sample**  
**Package insert & User Guide for the MEDICAL LABORATORY**

See the glossary of symbols at the end of the package insert

**INTENDED USE**

The vaginal FLOQSwabs® device for self-collection is a dry, sterile swab contained in a tube intended to be used for self-collection of vaginal samples for the detection of microorganisms (bacteria, fungi and viruses) responsible for sexually transmitted infections (STI). Once the vaginal sample has been collected, the tube allows transporting it from the collection site to the medical laboratory. In the laboratory, the samples can be analysed using diagnostic procedures for the detection of DNA of pathogenic agents.

**SUMMARY AND PRINCIPLES**

Sexually transmitted infections may be caused by different types of bacteria, fungi and viruses. Faced with such an extensive range of pathogens, indications by the physician and the reference guidelines can help the laboratory establish which analyses are appropriate to identify the aetiological agents of these infections (1). One of the routine procedures used in diagnosis of STI in women includes taking vaginal samples with a swab. This procedure can be carried out using the vaginal FLOQSwabs® device for self-collection.

The vaginal FLOQSwabs® device is a sterile flocked swab contained in a tube which allows collecting the vaginal sample both at a surgery and at home, as it has also been developed for self-collection. Self-collected vaginal swabs analysed using nucleic acid gene amplification methods are equivalent to clinician-collected swabs and, in some cases, more sensitive in detecting the pathogen responsible for the infection than the conventional bacterial culture methods (2,3), as the search for the nucleic acid does not require keeping the pathogen intact or vital. In particular, self-collected vaginal samples have provided better results than those collected by a medical professional in the search for *Chlamydia trachomatis* (CT) and *Neisseria gonorrhoeae* (NG) by means of molecular amplification (4).

The vaginal FLOQSwabs® device is composed of a plastic tube and a shaft with a flocked tip for collecting the sample. The shaft has a red line as reference marker that indicates where to position the fingers during self-collection. After collection, the sample must immediately be inserted in the transport tube and sent to the laboratory to perform the analysis.

**PRODUCT DESCRIPTION**

The vaginal FLOQSwabs® device for self-collection is composed of a plastic shaft for sample collection and a plastic tube for sample transport. One end of the shaft has a tip covered with Nylon fibres that act as a soft brush facilitating and improving collection of the biological material and the opposite end has a pink cap for closing the tube. The shaft has a red reference line that indicates where to position the fingers during the collection procedure for easy and correct gripping and handling of the device. Once the vaginal sample has been collected, the tube allows transporting it from the collection site to the analysis laboratory (Figure 1).

Figure 1 - Vaginal FLOQSwabs® device for self-collection



Figure 2 - Removable label

The tube has a removable barcode label with unambiguous numbering identical to that on the tube label. The removable barcode label can be detached and stuck onto the documentation provided to the analysis laboratory to guarantee traceability of the sample (Figure 2).

**REQUIRED MATERIALS BUT NOT INCLUDED**

Materials suitable for extraction, amplification and analysis of the nucleic acid of the pathogen searched for. For these procedures, consult the reference laboratory manuals.

**PRODUCT STORAGE**

It must be stored in its original package at a temperature between 2 and 30°C until the time of use. Do not incubate or freeze prior to use. If incorrectly stored, its effectiveness will be compromised.

**LIMITATIONS**

1. The appropriate environment for use of the device (collection in a hospital, surgery or home environment) must be defined by the healthcare professional based on the diagnostic investigation.
2. If the sample collected contains traces of blood (the sample will have a red or brown colour), it should be discarded and not used for analysis.
3. Vaginal samples collected with the swab at home and transported dry (without using a preservation culture medium) should not be analysed for the detection of STI using culture plate methods, as dry transport does not guarantee optimal survival of the strains involved.
4. For self-collection of vaginal samples in under-age patients, it is recommended to refer to the local, national and federal regulations.
5. The collection and transport conditions, the times and volume of the samples collected for clinical investigation are significant variables for obtaining reliable results.

**WARNINGS AND PRECAUTIONS FOR USE**

1. This product is for single use only; reuse may cause a risk of cross-contamination and/or inaccurate results.
2. A vaginal sample must not be self-collected during the menstrual cycle, pregnancy or in the three months after giving birth, and if the patient has recently had a gynaecological operation, or in case of unusual bleeding or pelvic pains.
3. In the 3 days prior to self-collection of a vaginal sample, use of vaginal ovules, creams or washes, vaginal contraceptives or condoms should be avoided.
4. In the 2 days prior to self-collection of a vaginal sample, sexual intercourse and ultrasound scans or gynaecological examinations should be abstained from.
5. Do not sterilize the device before use.
6. Do not repack for other commercial uses.
7. Not suitable for applications other than those for which it is intended.
8. Do not use if the device is visibly damaged (if the swab tip or shaft are broken).
9. The analysis method used for diagnostic purposes starting from vaginal samples collected with the vaginal FLOQSwabs® device for self-collection must be validated by the laboratory in advance.
10. In the laboratory, the use of this product in combination with preservation culture mediums, rapid diagnostic kits or diagnostic instruments must be validated by the user in advance.
11. All samples taken should be considered potentially infected, therefore, maximum caution is recommended when handling the sample.
12. Once the sample has arrived at the laboratory, wear protective gloves and other protective equipment suitable for handling clinical samples. Observe the biosafety recommendations when handling or analysing clinical samples.
13. After use, the swab and the transport tube must be disposed of in compliance with the laboratory regulations for infected waste.

**PRODUCT DETERIORATION**

Do not use the vaginal FLOQSwabs® device if: 1. the product shows visible signs of damage or contamination; 2. the product shows visible signs of sample leakage; 3. the product has passed its expiry date; 4. the tube is open; 5. or the product shows any other signs of deterioration.

**INSTRUCTIONS FOR USE****Self-collected vaginal samples**

For detailed instructions on self-collection of vaginal samples, refer to the **INSTRUCTIONS FOR USE FOR THE PATIENT** provided with the device.

**Vaginal sample transport and storage**

The vaginal sample collected with the vaginal FLOQSwabs® device can be transported in the swab tube. It is not necessary to dry the sample before inserting the swab in the tube.

For the detection of CT and NG by means of DNA amplification, the sample collected can be stored in the tube at ambient temperature (20-25°C) or refrigerated temperature (4°C) for up to 4 weeks.

For sample transport and handling, always comply with the national and federal regulations for transporting samples containing possible aetiological agents.

**Use in the laboratory for molecular analysis**

The vaginal sample collected with vaginal FLOQSwabs® device should be analysed as soon as it arrives at the laboratory. In any case, refer to the sample storage conditions included in the instructions for use of the molecular diagnostic kit used for the microorganism under examination.

The vaginal sample collected with vaginal FLOQSwabs® device must be diluted before amplification in an appropriate medium for molecular investigation. The swab must be diluted in a test tube containing PBS solution or another medium. The type and volume of medium will depend on the protocol used for the diagnostic molecular examination and the procedure validated by the laboratory.

**QUALITY CONTROL**

The vaginal FLOQSwabs® device for self-collection is subjected to various quality control tests before the product is released, including the test for the absence of DNase and RNase.

**PERFORMANCE CHARACTERISTICS**

The product has been validated on GeneXpert® molecular platform with Cepheid's Xpert® CT/NG kit, inoculating on the swab a suspension 10X LoD and 1X LoD of *Neisseria gonorrhoeae* (ATCC® 43069) and *Chlamydia trachomatis* (ATCC® VR822), and verifying the positivity for these two strains on swabs stored in their tube for 1, 2, 3 and 4 weeks at +4°C and +30°C.

To this end, after inoculation, the swab was diluted for 20 seconds in a test tube containing 2.5ml of sterile PBS of which 1ml was subsequently processed with the Cepheid's Xpert® CT/NG kit according to the instructions for use of the kit.

## RESULTS

The results obtained depend on appropriate and adequate collection of the sample, as well as on prompt transport and laboratory analysis. Incorrect analysis results may be obtained following collection errors, incorrect sample handling or storage, sample exchange, or presence in the sample of a number of organisms less than the analytical sensitivity threshold of the diagnostic test used.

In order to prevent incorrect results, strictly follow the instructions described in this leaflet and those provided by the manufacturer of the diagnostic kit used.

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UA.TR.116

#### BIBLIOGRAFIA

1. Magnus Unemo. Laboratory diagnosis of sexually transmitted infections, including human immunodeficiency virus. 2013. World Health Organization. Ronald Ballard, Catherine Ison, David Lewis, Francis Ndowa, Rosanna Peeling.
2. Catherine M W Stewart - Assessment of self taken swabs versus clinician taken swab cultures for diagnosing gonorrhoea in women: single centre, diagnostic accuracy study -BMJ 2012;345:e8107 doi: 10.1136/bmj.e8107 (Published 12 December 2012).
3. Masek BJ et al. *Performance of three nucleic acid amplification tests for detection of Chlamydia trachomatis and Neisseria gonorrhoeae by use of selfcollected vaginal swabs obtained via an Internet-based screening program.* Journal of Clinical Microbiology, 2009, 47(6):1663–1667.
4. Hobbs MM et al. *From the NIH: proceedings of a workshop on the importance of self-obtained vaginal specimens for detection of sexually transmitted infections.* Sexually Transmitted Diseases, 2008, 35(1):8–13.

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FLOQSwabs® are Copan Italia trademarks.  
 FLOQSwabs® are PATENTED WORLDWIDE No. PCT/EP2004/003392 and in the following countries: European Patent #1608268, US Patent #US8114027B2, Canadian Patent #2515205, Japanese Patent #2007-523663, Australian Patent #2004226798, New Zealand Patent #541560, Chinese Patent #101103931.

Copan Italia S.p.A.  
Via F. Perotti, 10  
25125 Brescia Italy  
Tel +39 030 2687211  
Fax +39 030 2687250

Email: [info@copangroup.com](mailto:info@copangroup.com)  
Website: [www.copangroup.com](http://www.copangroup.com)

North American Distributor:  
Copan Diagnostics Inc.  
26055 Jefferson Avenue  
Murrieta, CA 92562 USA  
Tel: 951-696-6957  
Fax: 951-600-1832

E-mail: [customerservice@copanusa.net](mailto:customerservice@copanusa.net)  
Website: [www.copanusa.com](http://www.copanusa.com)

Innovating Together™