Copan

HPC146R05 Date 2022.03

self
FLOQSwabs[®]

Instructions for Use

CE 0123

Self-Vaginal FLOQSwabs® with breakpoint
Instructions for use for the MEDICAL LABORATORY

INTENDED USE

FLOQSwabs® for self-collection of a vaginal sample is a dry, sterile swab contained in a tube intended to be used for self-collection of vaginal samples to test for microorganisms (bacteria, fungi and viruses) responsible for sexually transmitted infections (STI). Once the vaginal sample has been collected, the tube allows transporting it from the collection site to the medical laboratory. In the laboratory, the samples can be analysed using diagnostic procedures for the detection of DNA of pathogenic agents.

SUMMARY AND PRINCIPLES

Sexually transmitted infections may be caused by different types of bacteria, fungi and viruses. Faced with such an extensive range of pathogens, indications by the physician and the reference guidelines can help the laboratory establish which analyses are appropriate to identify the aetiological agents of these infections⁽¹⁾. One of the routine procedures used in diagnosis of STI in women includes taking vaginal samples with a swab. This procedure can be carried out using the vaginal FLOQSwabs® device for self-collection.

Self-vaginal FLOQSwabs® is a sterile flocked swab contained in a tube which allows vaginal self-collection in clinical setting. Self-collected vaginal swabs analysed using DNA amplification methods are equivalent to clinician-collected swabs and, in some cases, more sensitive in detecting the pathogen responsible for the infection than the conventional bacterial culture methods^(2,3), as the search for the nucleic acid does not require keeping the pathogen intact or vital. In particular, self-collected vaginal samples have provided better results than those collected by a medical professional in the search for *Chlamydia trachomatis* (CT), *Neisseria gonorrhoeae* (NG) and Human Papilloma Virus (HPV) by means of molecular amplification^(4,5). Using the vaginal sample self-collection procedure to search for HPV has moreover allowed optimizing the screening analyses for this pathogen type, considered closely related to the development of some forms or tumour of the uterine cervix⁽⁶⁾. The Self-vaginal FLOQSwabs® is composed of a plastic tube and a shaft with a flocked tip for collecting the sample. The shaft has a premoulded breakpoint with a red line. After collection, the sample must immediately be inserted in the transport tube and sent to the laboratory to perform the analysis.

PRODUCT DESCRIPTION

The Self-vaginal FLOQSwabs® with breakpoint is composed of a plastic shaft for sample collection and a plastic tube for sample transport. One end of the shaft has a tip covered with Nylon fibres and the opposite end has a cap that allows inserting the shaft in the tube and closing it. The shaft has a premoulded breakpoint at about 80mm or 60 mm from the tip based on the catalogue number. The breaking point allows the break of the swab inside a test tube if required by the laboratory procedure. The breakpoint is identified with a red mark (Figure 1). This red mark also indicates where to position the fingers during the collection procedure for easy and correct gripping and handling of the device. Once the vaginal sample has been collected, the tube allows transporting it from the collection site to the analysis laboratory.

Some catalogue numbers have a label with a removable barcode on the tube. The removable barcode has numbering identical to the one on the tube label. The removable barcode label can be detached and stuck onto the documentation provided to the analysis laboratory to guarantee traceability of the sample (Figure 2).

Figure 1 - FLOQSwabs® with breakpoint for self-collection of a vaginal sample



Figure 2 - Removable label

FLOQSwabs® for self-collection of a vaginal sample is available in different product configuration described in the table below:

Catalogue Number	Description	Pack size
552C.80	One regular size applicator swab with 80mm breaking point, in tube with red cap	100 units per box/ 6 boxes x 100 units in a cardboard
552C.80PB	One regular size applicator swab with 80mm breaking point, in tube with red cap and removable barcode	100 units per box/ 6 boxes x 100 units in a cardboard
5E089N	One regular size applicator swab with 60mm breaking point, in tube with white cap and removable barcode	100 units per box/ 6 boxes x 100 units in a cardboard

REQUIRED MATERIALS BUT NOT INCLUDED

Materials suitable for extraction, amplification and analysis of the DNA of the pathogen searched for. For these procedures, consult the reference laboratory manuals.

PRODUCT STORAGE

It must be stored in its original package at a temperature between 2 and 30°C until the time of use. Do not incubate or freeze prior to use. If incorrectly stored, its effectiveness will be compromised.

LIMITATIONS

1. Vaginal samples collected and transported dry (without using a preservation culture medium) should not be analysed to search for STI using culture plate methods, as dry transport does not guarantee optimal survival of the strains involved.
2. For self-collection of vaginal samples in under-age patients, it is recommended to refer to the local, national and federal regulations.
3. The collection and transport conditions, the times and volume of the samples collected for clinical investigation are significant variables for obtaining reliable results.
4. The FLOQSwabs® vaginal swab has a premoulded breakpoint. If the laboratory procedure requires breaking it in a test tube, before proceeding with breaking, check that the height of the test tube is compatible with the height of the breakpoint (from the tip to the breakpoint). If test tube different from Copan 12x80 mm test tube are used and the length of the shaft (from the tip to the breakpoint) inside the test tube does not fit the height of the test tube, the breaking mechanism described in these instructions for use may not be applicable and it may be impossible to close the test tube.

WARNINGS AND PRECAUTIONS FOR USE

1. The appropriate environment for use of the device (collection in a hospital, surgery or home environment) must be defined by the healthcare professional based on the diagnostic investigation.
2. This product is for single use only; reuse may cause a risk of cross-contamination and/or inaccurate results.
3. The vaginal self-collection should not be performed: 1) during pregnancy; 2) in the three months after giving birth; 3) if the patient recently had a gynaecological operation; 4) in case of unusual bleeding or pelvic pains. Specimen collection during the menstrual cycle has not been evaluated.
4. In the 3 days prior to self-collection of a vaginal sample, use of vaginal ovules, creams or washes, vaginal contraceptives or condoms should be avoided.
5. In the 2 days prior to self-collection of a vaginal sample, sexual intercourse and ultrasound scans or gynaecological examinations should be abstained from.
6. Do not sterilize the device before use.
7. Do not repack for other commercial uses.
8. Not suitable for applications other than those for which it is intended.
9. Do not use if the device shows visible signs of sample contamination or leakage upon arrival at the laboratory.
10. The analysis method used for diagnostic or screening purposes starting from vaginal samples collected with the vaginal FLOQSwabs® device for self-collection must be validated by the laboratory in advance.
11. In the laboratory, the use of this product in combination with preservation culture mediums, rapid diagnostic kits or diagnostic instruments must be validated by the user in advance.
12. All samples taken should be considered potentially infected, therefore, maximum caution is recommended when handling the sample.
13. Once the sample has arrived at the laboratory, wear protective gloves and other protective equipment suitable for handling clinical samples. Observe the biosafety recommendations when handling or analysing clinical samples.
14. After use, the swab and the transport tube must be disposed of in compliance with the laboratory regulations for infected waste.

PRODUCT DETERIORATION

Before handling FLOQSwabs® to the patient for self-collection, check that the product is intact and properly preserved. Do not use if: 1. the product shows visible signs of damage (e.g. if the swab tip or shaft is broken) or contamination, 2. the product has passed its expiry date, 3. the tube is open (seal not intact), 4. or the product shows any other signs of deterioration.

INSTRUCTIONS FOR USE

Self-collection of vaginal samples

For detailed instructions on self-collection of vaginal samples, refer to the **INSTRUCTIONS FOR USE FOR THE PATIENT** provided with the device.

Vaginal sample transport and storage

The vaginal sample self-collected with FLOQSwabs® can be transported in the swab tube. It is not necessary to dry the sample before inserting the swab in the tube.

Please refer to the downstream diagnostic test kit for the specific transport and storage claims of the sample.

Use in the laboratory for molecular analysis

The vaginal sample collected with FLOQSwabs® should be analysed as soon as it arrives at the laboratory.

The vaginal sample collected with FLOQSwabs® must be diluted before amplification in an appropriate medium for molecular investigation. The type and volume of medium will depend on the protocol used for the diagnostic kit and the procedure validated by the laboratory.


Procedure to break the FLOQSwabs® shaft in a test tube (Fig. 3)

If the laboratory procedure requires breaking the swab inside a test tube, break it at the premoulded breakpoint as follows:

NOTES: The instructions below describe the operating procedure to break FLOQSwabs® inside a standard Copan 12x80mm test tube. The breaking mechanism with other test tubes and/or swabs different from those indicated must be evaluated by the user in advance and may require different or additional operating instructions and precautions from those described below.

- Hold the test tube in one hand with the test tube opening facing away from your face.
- Grip the end of the shaft with the thumb and forefinger of the other hand.
- Rest the part of the shaft with the breakpoint against the edge of the test tube.
- Bend the swab shaft at an angle of 180° so that it breaks in correspondence to the breakpoint. If necessary, gently turn the swab shaft until it is completely broken off and remove the upper part of the shaft.
- Discard the broken off part of the swab shaft.
- Position the cap on the test tube and close it.

Fig. 3 - Procedure to break the self-vaginal FLOQSwabs® shaft in a Copan 12x80mm test tube

QUALITY CONTROL

The self-vaginal FLOQSwabs® is subjected to various quality control tests before the product is released, including the test for the absence of DNase and RNase.

PERFORMANCE CHARACTERISTICS

For the CT/NG test, the product was validated on GeneXpert® molecular platform with Xpert® CT/NG kit (Cepheid), inoculating on the swab a suspension 10X LoD and 1X LoD of *Neisseria gonorrhoeae* (ATCC® 43069) and *Chlamydia trachomatis* (ATCC® VR822). After inoculation, the swab was broken inside the test tube of the Cepheid's Xpert® CT/NG kit and the sample processed according to the instructions for use of the kit. CT/NG DNA stability onto dry swab was tested with Cepheid's Xpert® CT/NG kit at time zero and after 4 weeks stored at room temperature (20-25°C) and refrigerated temperature (4°C) by eluting the swabs in 2.5 ml of phosphate- buffered solution.

For the HPV test, the product was validated on cobas® 4800 molecular platform with cobas® HPV test kit (Roche Molecular Systems, Inc). The swab was immersed in a suspension of HPV18 positive cells for 10 seconds at two different concentrations: 10X LoD (400 cell/mL) and 1X LoD (40 cell/mL). After inoculation, the swab was placed in 5ml of PreservCy® Solution (Hologic, Inc.) and swirled in the solution for 20 seconds. The sample was then processed with Roche cobas® HPV test kit according to the instructions for use of the kit. The same procedure was followed to test the stability of the sample at time zero and after two weeks stored at -20°C and +50°C by inoculating the swab with a suspension of HPV16 positive cells with a concentration of 300 cells/mL.

The product code 5E089N, was validated using clinical specimens tested with:

- 1) BD Viper LT and BD COR™ automated platforms with the BD Onclarity HPV Assay (Beckton Dickinson). After the specimen collection, in the laboratory the self-collected swabs were transferred to BD HPV LBC Self-Collection diluent (3 mL) tubes and processed with BD Onclarity HPV Assay on both BD Viper LT and BD COR™ automated platforms according to the manufacturer instructions.
- 2) Ampliquality HPV Type Express (AB Analytica srl). After the specimen collection, the swab was immersed in 200 µl of the NucleoSpin 8 Tissue (Macherey-Nagel) lysis buffer for 2 hours at 56°C to extract the viral DNA as per the kit manufacturer's instructions. The sample was then processed with test kit Ampliquality HPV Type Express (AB Analytica) according to the instructions for use of the kit.

RESULTS

The performance results or the stability of the specimen obtained depend on appropriate and adequate collection of the sample, prompt transport and laboratory analysis as well as the downstream assay in use. Incorrect analysis results may be obtained following collection errors, incorrect sample handling or storage, or presence in the sample of a number of microorganisms less than the analytical sensitivity threshold of the diagnostic test used.

In order to prevent incorrect results, strictly follow the instructions described in this leaflet and those provided by the manufacturer of the diagnostic kit used.

Please refer to symbol table at the end of the instructions for use

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UA.TR.116

BIBLIOGRAPHY

1. Magnus Unemo. Laboratory diagnosis of sexually transmitted infections, including human immunodeficiency virus. 2013. World Health Organization. Ronald Ballard, Catherine Ison, David Lewis, Francis Ndowa, Rosanna Peeling.
2. Catherine M W Stewart - Assessment of self taken swabs versus clinician taken swab cultures for diagnosing gonorrhoea in women: single centre, diagnostic accuracy study -BMJ 2012;345:e8107 doi: 10.1136/bmj.e8107 (Published 12 December 2012)
3. Masek BJ et al. *Performance of three nucleic acid amplification tests for detection of Chlamydia trachomatis and Neisseria gonorrhoeae by use of selfcollected vaginal swabs obtained via an Internet-based screening program.* Journal of Clinical Microbiology, 2009, 47(6):1663–1667.
4. Hobbs MM et al. *From the NIH: proceedings of a workshop on the importance of self-obtained vaginal specimens for detection of sexually transmitted infections.* Sexually Transmitted Diseases, 2008, 35(1):8–13.
5. Perderson HM et al (2018) *Implementation considerations using HPV self-collection to reach women under-screened for cervical cancer in high-income settings.* Curr Oncol. 25(1):e4-e7
6. Arbyn M et al (2018) *Detecting cervical precancer and reaching underscreened women by using HPV testing on self-samples: updated meta-analyses.* BMJ. 363: k4823.

Index of Symbols

Symbol / Simbolo / Symbol / Symbole / Símbolo / Símbolos / Symbol / Symbol / Символ	Meaning / Significato / Bedeutung / Sens / Signification / Significado / Betekenis / Betydelse / Значення
	Manufacturer / Fabbicante / Hersteller / Fabricant / Fabricante / Fabricante / Fabrikant / Tilverkare / Виробник
CE 0123	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificación del organismo notificado / Identificação do organismo notificado / Identificatienummer aangemelde instantie / Det anmeldte organets identifikationskod / Ідентифікаційний номер уповноваженого органу
	Sterilized using ethylene oxide / Sterilizzato usando ossido di etilene / Sterilisiert mit Äthylenoxid / Stérilisé à l'aide d'oxyde d'éthylène / Esterilizado por óxido de etileno / Esterilizado utilizando óxido de etileno / Gesteriliseerd met ethyleenoxide / Steriliserat med etylenoxid / Стерилізований завдяки етиленоксиду
	Do not reuse / Non riutilizzare / Nicht zur Wiederverwendung / Ne pas réutiliser / No reutilizar / Não voltar a utilizar / Niet opnieuw gebruiken / Ateranvänd inte / Не використовувати повторно
	Catalogue number / Numero di catalogo / Bestellnummer / Référence du catalogue / Número de catálogo / Número do catálogo / Catalogusnummer / Katalognummer / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Limites de temperatura / Temperaturgrenzen / Temperaturgränser / Температурні обмеження
	Use by / Utilizzare entro / Verwendbar bis / Utiliser jusque / Fecha de caducidad / Prazo de validade / Te gebruiken tot / Sista förbrukningsdag / Використати до
	Consult Instructions for Use / Consultare le istruzioni per l'uso / Gebrauchsanweisung beachten / Consulter les instructions d'utilisation / Consulte las instrucciones de uso / Consultar as instruções de utilização / Raadpleeg de gebruiksaanwijzing / Läs bruksanvisningen / Користуватись інструкціями по використанню
	Batch code (Lot) / Codice del lotto (partita) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código de lote (Lote) / Código do lote (Lote) / Partijcode (Lotnummer) / Partikod (Lot) / Код серії (партії)
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> test / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contenido suficiente para <n> pruebas / Contém o suficiente para <n> testes / Bevat voldoende voor <n> testen / Innehåller tillräckligt för <n> tester / Придатний для проведення <n> кількості випробувань
	Do not use if package is damaged / Non utilizzare in caso di confezionamento danneggiato / Bei Beschädigung der Verpackung nicht verwenden / Ne pas utiliser si l'emballage est abîmé / No utilizar en caso de paquete dañado / Não utilizar se a embalagem estiver danificada / Niet gebruiken wanneer de verpakking beschadigd is / Använd inte vid skadad förpackning / Не використовувати, якщо упаковка пошкоджена

FLOQSwabs® are Copan Italia trademarks.
 FLOQSwabs® are PATENTED WORLDWIDE No. PCT/EP2004/003392 and in the following countries: European Patent #1608268, US Patent #US8114027B2, Canadian Patent #2515205, Japanese Patent #2007-523663, Australian Patent #2004226798, New Zealand Patent #541560, Chinese Patent #101103931.

Copan

Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia, Italy
Tel +39 030 2687211
Fax +39 030 2687250

Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:

Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com