**AGITADOR MAGNÉTICO
DIGITAL MAGNETIC STIRRER
AGITATEUR MAGNETIQUE NUMERIQUE**

Serie / Series / Séries / RSLAB-2NC, RSLAB-2C

RS Lab

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This manual should be available for all users of these equipments. To get the best results and a higher duration of this equipment it is advisable to read carefully this manual and follow the processes of use.

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
INSTRUCTIONS ON ENVIRONMENT PROTECTION

At the end of its life cycle, please, does not dispose of this equipment by throwing it in the usual garbage; hand it over a collection point for the recycling of electrical and electronic appliances. It does not contain dangerous or toxic products for humans but a non-adequate disposal would damage the environment.

The materials are recyclable as mentioned in its marking. By recycling material or by other forms of re-utilization of old appliances, you are making an important contribution to protect our environment.

Please inquire at the community administration for the authorized disposal location.

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INDEX OF CONTAINS

1. Safety instructions	15	b. Residual heat warning	20
2. Proper use	16	c. Setting the safety temperature	20
3. Inspection	16	9. Function: stirring	
4. Trial run	16	a. Basic stirring	21
5. Control and display		b. Stirring bar breakaway monitoring	21
a. Control elements	17	10. Remote control	21
b. Display	18	11. Faults	22
6. Initialization	18	12. Maintenance and cleaning	22
7. Operating modes		13. Associated standards and regulations	23
a. Operating mode A	19	14. Technical data	23
b. Operating mode B	19	15. Main spare parts	24
c. Operating mode C	19	16. Accessories	25
8. Function: heating (RESLAB-2C)			
a. Working with external temp. sensor	20		

Welcome to the Digital magnetic stirrer RSLAB-2NC (no heating) and RSLAB-2C (heating) user manual. Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

When help needed, you can always contact the service department of manufacturer or your supplier for technical support. Please provide the customer care representative with the following information: Serial number, description of problem, methods and procedures adopted to resolve the problems and your contact information.

RSLA stirrers are warranted to be free from defects in materials and workmanship under normal use and service for a period of 12 months from the date of invoice. The warranty is extended only to the original purchaser. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claims under the warranty please contact your local dealer. You may also send the instrument direct to our works, enclosing the invoice copy and giving reasons for the claim.

1- SAFETY INSTRUCTIONS**WARNING!**

Read the operating instructions carefully before use the instrument.

Ensure that only trained staff work with the instrument.

Forbid to heat the substances with low-burning point or easy-volatile.

Risk of burn!


- ◆ Caution when touching the housing parts and the heating plate. The heating plate can reach temperatures of 340 °C.
- ◆ Pay attention to the residual heat after switching off.

Protective ground contact!

- ◆ Make sure that socket is earthed (protective ground contact) before use.
- ◆ When work, wear the personal guard to avoid the risk from splashing and evaporation of liquids and release of toxic or combustible gases.
- ◆ Set up the instrument in a spacious area on a stable, clean, non-slip, dry and fireproof surface, do not operate the instrument in explosive atmospheres, with hazardous substances or under water.
- ◆ Gradually increase the speed, reduce the speed if the stirring bar breakaway because of too high speed and if the instrument is not running smoothly, or if the container moves on the base plate.
- ◆ Temperature must always be set to at least 25 °C lower than the fire point of the media used.
- ◆ Beware of hazards due to flammable materials or media with a low boiling temperature, overfilling of media and unsafe container.
- ◆ Process pathogenic materials only in closed vessels.
- ◆ Check the instrument and accessories before handle it for damaged component. Safe operation is only guaranteed with the accessories described in the "Accessories" chapter. Accessories must be securely attached to the device and cannot come off by themselves. Always disconnect the plug before fitting accessories.
- ◆ Ensure that the external temperature sensor is inserted in the media to a depth of at least 20 mm.
- ◆ When using metal vessels, do not place the temperature sensors on the bottom of the vessel. Placing sensors on the vessel bottom can cause excessively high temperatures to be measured especially in media which have poor conductivity. The tip of the measuring sensor must be at least 5 mm from the vessel bottom; a distance of 10 mm is ideal.
- ◆ The instrument can only be disconnected from the main power supply by pulling out the main plug or the connector plug.
- ◆ The voltage started on the label must correspond to the main power supply.
- ◆ Ensure that the main power supply cable does not touch the heating base plate. Do no

15- MAIN SPARE PARTS**N° Description**

- 1 Work plate
- 2 Hole for external temp. sensor holder
- 3 Mains switch
- 4 Housing
- 5 Ventilation
- 6 Control panel
- 7 External temperature sensor's socket
- 8 Remote control socket
- 9 Power supply socket
- 10 Main circuit board
- 11 Driver circuit board
- 12 Motor
- 13 Turntable
- 14 Internal temperature sensors
- 15 Permanent magnet
- 16 Heating parts
- 17 Protective plate
- 18 Rear panel
- 19 Robber bush
- 20 Interface circuit board

- ◆ Before using other method for cleaning or decontamination, the user must ascertain with the manufacturer that this method does not destroy the instrument.
- ◆ Unit should be clean and fitted in the original packing or a suitable packing for transport.
- ◆ The enamel makes the heating plate easier to care for and more resistant to acids and bases. Because of it, however, the heating plate is also more susceptible to extreme fluctuations in temperature and the force of impact. This can result in cracks forming or the coating flaking off.

13- ASSOCIATE STANDARDS AND REGULATIONS

Construction in accordance with the following safety standards:

EN 61010-1
 UL 3101-1
 CAN/CSA C22.2 (1010-1)
 EN 61010-2-10

Construction in accordance with the following EMC standards:

EN 61326-1

Associated EU guidelines:

EMC guidelines: 89/336/EWG
 Instrument guidelines: 73/023/EWG


14- TECHNICAL DATA

Voltage: 200-240/100-120VAC
Frequency: 50/60 Hz
Power: 550W (con calefacción) / 50W
Stirring point position quantity: 1
Max. stirring quantity (H₂O): 20 Litres
Max. magnetic bar: 80mm
Motor type: Brushless motor with external rotor
Max. power input of motor: 18W
Max. power output of motor: 10W
Speed range: 100-1500rpm
Rotary speed display: LCD
Speed display accuracy: 1rpm
Ø Of the heating plate: 135mm
Heating power: 500W
Heating plate material: Stainless steel Porcelain enamel
Control accuracy of heating temperature with temperature sensor: 0,2°C

Heating rate (1 l water): 6Kw/min
Temperature range: RT-340°C
Temperature display: LCD
Temperature display accuracy: 0,1°C
Control accuracy of heating temperature: 1°C
The safety temp. rage of the hotplate: 100-350°C
Temperature sensor: PT1000

Dimensions: 280x160x85mm
Weight: 2,8kg
Permitted ambient temperature: 5-40°C
Permitted relative humidity: 80%
Interface RS232: Si
Protection class acc. to DIN 60529: IP42

- cover the device.
- ◆ The instrument may only be opened by experts.
- ◆ Keep away from high magnetic field.
- ◆ Observe the minimum distances between the devices, between the instruments and the wall and above the assembly (min. 100 mm).

2- PROPER USE

The instrument is designed for mixing and/or heating liquids in schools, laboratories or factories. This device is not suitable for using residential areas or other constraints mentioned in chapter 1.

3- INSPECTION

Unpack the equipment carefully and check for any damages which may have arisen during transport. If it happens, please contact manufacturer for technical support.

If there is any apparent damage to the system, please do not plug it into the power line.

The RSLAB-2NC y RSLAB-2C includes the following items:


Main unit: 1 unit
 Power cable: 1 unit
 Stirrer bar: 1 unit
 User manual: 1 unit

4- TRIAL RUN

- ◆ Make sure the required operating voltage and power supply voltage match.
- ◆ Ensure the socket must be earthed reliably.
- ◆ Ensure the power be off.
- ◆ Plug in the power cable, ensure the power be on and begin initializing.
- ◆ Add the medium into the vessel with a stirring bar.

- ◆ Put the vessel on the workplate.
- ◆ Set the rated stirring speed and start stirring.
- ◆ Observe the stirring bar and LCD display.
- ◆ Set the rated temperature and start heating.
- ◆ Observe the real temperature on LCD display.
- ◆ Stop the heating and stirring functions.
- ◆ If these operations above are normal, the device is ready to operate. If these operations are not normal, the device may be damaged during transportation, please contact manufacturer for technical support.

5- CONTROL AND DISPLAY

1. LED display
2. Mode key
3. Function key
4. Heating Knob (only 2C model)
5. Stirring Knob
6. LCD display
7. On/Off

Items	Descriptions
Stirring Knob	Set the rated rotary speed. The function "stirring" is switched On or OFF via push the knob
Heating Knob (RSLAB-2C)	Set the rated temperature. The function "heating" is switched ON or OFF via push the knob
Mode Key	Shift among the three operating modes: A, B and C mode
Function Key	Setting all parameters of the device
LCD	LCD displays the real working state and all settings
LED Heat (RSLAB-2C)	When the heating function is switched On, the LEAD heat is lit
LED Stir	When the stirring function is switched On, the LEAD stir is lit
Power Switch	Switch On or OFF the instrument

11- FAULTS

1. **Device can't be power ON:**
 - Check whether the power cable is plugged
 - Check whether the fuse is broken or loose
2. **Fault in power on self-test:**
 - Switch OFF the unit, then switch on and reset the instruments to factory default setting
3. **Temperature cannot reach set point:**
 - Check whether the safety temperature value is set too low
4. **Stir speed cannot reach the set point:**
 - Check whether the stirring bar breakaway monitoring is ON
 - Excessive medium viscosity may cause abnormal speed reduction of the motor
5. **Heating or stirring cannot be started when push the temperature control knob:**
 - Check whether the unit is in the mode setting program
 - Exit from the mode setting and restart the heating/stirring functions
6. **Unit cannot be powered off when switched off:**
 - Check if the residual heat warning functions is still ON and hotplate temperature is above 50 °C (The LCD still work and "Hot" flash)
 - Turn off residual heat warning function before powered OFF the unit

If these faults are not resolved, please contact manufacturer for technical support.

12- MAINTENANCE AND CLEANING

- ◆ Proper maintenance can keep instruments working in a good state and lengthen its life-time.
- ◆ Be careful not spray the cleaners into the instrument when cleaning.
- ◆ Unplug the power line when cleaning.
- ◆ Only use cleaners that we advise as below:
 - Dyes: isopropyl alcohol
 - Construction materials: water containing tenside / Isopropyl alcohol
 - Cosmetics: water containing tenside / Isopropyl alcohol
 - Foodstuffs: water containing tenside
 - Fuels: water containing tenside
- ◆ Wear proper protective gloves during cleaning of the instrument.

9- FUNCTION: STIRRING

a. Basic stirring


- ◆ The function “stirring” is switched on or off via the rotary stirring knob. The motor speed is set on the knob (100 to 1500 rpm in steps of 10 rpm).
- ◆ When both functions heating and stirring are switched on and those above operations are done, the LCD will shift to the speed value and come back to the temperature value in about 5 seconds.

b. Stirring bar breakaway monitoring

- ◆ Breakaway can occur when stirring at high speeds, when large stirring bars are being driven at speeds close to the maximum magnetically torque or sample viscosity changes.
- ◆ Stirring bar breakaway monitoring function can be set at the rotary speed range between 300 and 1500 rpm.
- ◆ In the A or B mode, the stirring bar breakaway monitoring can be set, as follows:
- ◆ Press continuously the function setting key till the LCD displays “br”.
- ◆ Press the function setting key to switch on or off the stirring bar breakaway monitoring function, and the LCD displays “On” or “OFF”.
- ◆ Press continuously the function setting key till “set” disappear.
- ◆ When stirring bar breakaway is detected, stirring is stopped immediately to enable the stirring bar to be picked up. After this, stirring starts again and rotary speed increases up to the set speed. If a further stirring bar breakaway occurs within the next 3 minutes of stirring at a constant set speed, the stirring bar will be picked up again and stirring will restart but this time at the set speed reduced by 100 rpm. The speed will be reduced each time in 100 rpm steps until the minimum of 200 rpm is reached. In this case you should select new test conditions.


IMPORTANTE: You must calibrate and test the stirring bar breakaway monitoring function for your particular configuration taking into account the stirring bars, speed, vessel and medium you are using.

10- REMOTE CONTROL

Connection PT 1000 RS 232C Socket

The device can be controlled from an external PC (using the dedicated software) via the RS 232 C serial interface fitted to the unit. Data communication from laboratory instrument to computer is only possible on demand of the computer.

Characters	Descripción
Temp et C°	Display temperature when the heating function is switched ON
Stir	Display stirring state when the stirring function is switched ON
Mode	Display the setting mode in order of A, B and C after pressing mode key
Set	Press setting key to set function, and display Set characters
Hot	Display hot warning if the heating plate temperature is above 50 °C after switching OFF the heating function
Probe	Display when using external probe
Remote	Display in case of remote control
Err	Display in case of error happening
Rated value	Display value in case of heating or stirring function switching ON
Real value	

If both heating and stirring functions have been started at the same time, the display of function “heating” always has higher priority. If in this case speed is changed via the stirring knob, it displays stirring speed and reverses to temperature in the duration on 5 seconds

6- INITIALIZATION

- ◆ Turn on instrument by the power switch.
- ◆ The “heat” and “stir” LEDs and the LCD display flash about 3 times.
- ◆ LCD displays the workplate safe temperature limit in “SAFE xxx °C”.
- ◆ LCD displays residual heat warning “rES On/OFF”
- ◆ LCD displays stirring bar breakaway monitoring “br On/OFF”
- ◆ LCD displays operating modes (A, B and C)
- ◆ LCD displays the “Probe” characters in case of connecting with external temperature sensors.

◆ RSLAB-2NC and RSLAB-2C are supplied with the following default settings:

- Operating mode: A
- Setting temperature: 25 °C
- SAFE: 350 °C
- Setting rotary speed: 100 rpm
- Residual heat warning: On
- Stirring bar breakaway monitoring: OFF

◆ Reset the parameters to the default settings, as follows:

Turn off the power switch, have the keys “Set” and “Mode” pressed down at the same time and switch on the power switch. Release buttons after 5 seconds. The unit is now set to factory default settings.

7- OPERATING MODES

The stirrer RSLAB-2NC y RSLAB-2C can be operated in the following three operation modes:

a. Operation mode A

- ◆ Heating and stirring functions.
- ◆ Users can control on the front panel or by remote PC.
- ◆ Residual temperature warning, safe temperature limit and stirring bar breakaway monitoring functions

b. Operation mode B

- ◆ Work with external temperature sensor.
- ◆ Heating and stirring functions.
- ◆ Users can control on the front panel or by remote PC.
- ◆ Residual temperature warning, safe temperature limit and stirring bar breakaway monitoring functions.
- ◆ Maintain the state of the heating and stirring functions before the last switching off.

c. Operation mode C

- ◆ Maintain the state of the heating and stirring functions before the last switching off.
- ◆ Workplate temperature, safe temperature limit and speed limit can't be set in this mode.
- ◆ Users can control on the front panel or by remote PC.
- ◆ Residual temperature warning, safe temperature limit and stirring bar breakaway monitoring functions.

Mode changing and parameter setting is forbidden when the instrument is running. Switch off “stir” and “heat” control knob in advance before change mode.

8- FUNCTION: HEATING (RSLAB-2C)

The stirrer hotplate is kept a constant temperature by a digital control circuit which has two separate safe circuits. The hotplate temperature can also be monitored from a separate adjustable safety circuit. The two temperature sensors (Pt 1000) internal for temperature control are built into the hotplate. The single external Pt 1000 can monitor the temperature of sample.

If this case, external temperature sensor must be plugged in before the device is switched on. Ensure the external temperature sensors connecting firmly before heating.

Set the temperature via turning the temperature control knob slowly to the ideal scale. When heating function is switched on, the LED “heat” will light up and the LCD will display the real temperature.

The rated temperature will be displayed on the right-hand side of the LCD as well as Temp and °C characters.

The heating function is switched on or off via pushing the heating knob.

a. Working with external temperature sensor.

In the B mode, if the external temperature sensor is not plugged in, “Probe” will flash. If plugged in, “Probe” will be shown on the digital display to indicate that the sensor is operating. The external temperature sensors' value and actual temperature are displayed.

Compared with the temperature control of the hotplate, the external temperature sensor can control the media temperature precisely.

b. Residual heat warning (HOT)

In order to prevent the risk of burns from a hotplate, RSLAB-2C has a residual heat warning function. When the heating function is switched off and the hotplate temperature is still above 50 °C, “hot” will flash to warn that there is a hazard of burns from the hotplate. When the heating plate temperature drops to below 50 °C, the unit will automatically switch off.

- ◆ In the A or B mode, users can set the residual heat warning function, as follow:
- ◆ Press continuously the function setting key till the LCD displays “rES”.
- ◆ Press the temperature regulation key to switch on or off the residual heat warning function, and LCD displays “On” or “Off”.
- ◆ Press continuously the Function setting key till “set” disappear.
- ◆ When the plug is pulled out, the residual heat warning function cannot be run.

c. Setting the safety temperature

- ◆ In the A or B mode, the safety temperature can be set, as follows:
- ◆ Press continuously the function setting key till the LCD displays “SAFE”.
- ◆ Adjust the temperature setting key to set the safety temperature from 100 to 350 °C.
- ◆ Press continuously the function setting key till “set” disappear.