
LANGE 

DOC022.98.90251

sensION™ + EC7

03/2013, Edition 3

User Manual
Bedienungsanleitung
Manuale dell'utente
Manuel d'utilisation
Manual del usuario
Manual do utilizador
Uživatelská příručka
Brugsanvisning
Gebruikershandleiding
Instrukcja obsługi
Bruksanvisning
Käyttäjän käsikirja
Ръководство на потребителя
Használati útmutató
Manual de utilizare
Naudotojo vadovas
Руководство пользователя
Kullanım Kılavuzu
Návod na obsluhu
Navodila za uporabo
Korisnički priručnik
Εγχειρίδιο χρήσης
Kasutusjuhend

English.....	3
Deutsch.....	14
Italiano.....	26
Français.....	38
Español.....	50
Português.....	62
Čeština.....	74
Dansk.....	85
Nederlands.....	96
Polski.....	108
Svenska.....	120
Suomi.....	131
български.....	142
Magyar.....	154
Română.....	165
lietuvių kalba.....	176
Русский.....	187
Türkçe.....	199
Slovenský jazyk.....	210
Slovenski.....	221
Hrvatski.....	232
Ελληνικά.....	244
eesti keel.....	256

Specifications

Specifications are subject to change without notice.

Specification	Details
Dimensions	35 x 20 x 11 cm (13.78 x 7.87 x 4.33 in.)
Weight	1100 g (2.43 lb)
Meter enclosure	IP42
Power requirements (external)	100–240 V, 0.4 A, 47-63 Hz
Meter protection class	Class II
Storage temperature	–15 to +65 °C (5 to +149 °F)
Operating temperature	0 to 40 °C (41 to 104 °F)
Operating humidity	< 80% (non-condensing)
Measuring error (± 1 digit)	Conductivity: ≤ 0.1%, Salinity: ≤ 0.5%, temperature: ≤ 0.2 °C (0.36 °F)
Reproducibility (± 1 digit)	Conductivity: ± 0.1%, Salinity: ± 0.1 %, temperature: ± 0.1 °C (0.18 °F)
Connections	Conductivity probe with built-in Pt1000 sensor: telephonic connector; magnetic stirrer: RCA connector
Temperature correction	Manual, Pt 1000 temperature probe (A.T.C.), NTC 10 kΩ probe
Measurement display lock	Continuous measurement, by stability
Display	Liquid crystal, backlit, 128 x 64 dots
Keyboard	PET with protective treatment
Certification	CE

General information

Revised editions are found on the manufacturer's website.

Safety information

NOTICE

The manufacturer is not responsible for any damages due to misapplication or misuse of this product including, without limitation, direct, incidental and consequential damages, and disclaims such damages to the full extent permitted under applicable law. The user is solely responsible to identify critical application risks and install appropriate mechanisms to protect processes during a possible equipment malfunction.

Please read this entire manual before unpacking, setting up or operating this equipment. Pay attention to all danger and caution statements. Failure to do so could result in serious injury to the operator or damage to the equipment.

Make sure that the protection provided by this equipment is not impaired. Do not use or install this equipment in any manner other than that specified in this manual.

Use of hazard information

▲ DANGER

Indicates a potentially or imminently hazardous situation which, if not avoided, will result in death or serious injury.

▲ WARNING

Indicates a potentially or imminently hazardous situation which, if not avoided, could result in death or serious injury.

▲ CAUTION

Indicates a potentially hazardous situation that may result in minor or moderate injury.

NOTICE

Indicates a situation which, if not avoided, may cause damage to the instrument. Information that requires special emphasis.

Precautionary labels

Read all labels and tags attached to the instrument. Personal injury or damage to the instrument could occur if not observed. A symbol, if noted on the instrument, will be included with a danger or caution statement in the manual.

	<p>This symbol, if noted on the instrument, references the instruction manual for operation and/or safety information.</p>
	<p>Electrical equipment marked with this symbol may not be disposed of in European public disposal systems after 12 August of 2005. In conformity with European local and national regulations (EU Directive 2002/98/EC), European electrical equipment users must now return old or end-of-life equipment to the Producer for disposal at no charge to the user.</p> <p>Note: For return for recycling, please contact the equipment producer or supplier for instructions on how to return end-of-life equipment, producer-supplied electrical accessories, and all auxiliary items for proper disposal.</p>

Product overview


The sensION™+ meters are used with probes to measure various parameters in water.

The sensION™+ EC7 meter measures conductivity, salinity and temperature.

Product components

Refer to [Figure 1](#) to make sure that all components have been received. If any items are missing or damaged, contact the manufacturer or a sales representative immediately.

Figure 1 Meter components





<p>1 Standard solutions (147 µS/cm, 1413 µS/cm and 12.88 mS/cm)</p>	<p>5 Rod with o-ring</p>
<p>2 Calibration beakers (with magnetic bar inside)</p>	<p>6 Probe (included with kits only)</p>
<p>3 Probe holder</p>	<p>7 Meter</p>
<p>4 Power supply</p>	

Installation

Assemble the probe holder

Follow the numbered steps to assemble the probe holder and to connect the magnetic stirrer.

Connect to AC power

⚠ DANGER


Electrocution hazard. If this equipment is used outdoors or in potentially wet locations, a Ground Fault Circuit Interrupt (GFCI/GFI) device must be used to connect the equipment to its main power source.

The meter can be powered by AC power with the universal power adapter.

1. Select the correct adapter plug for the power outlet from the adapter kit.
2. Connect the universal power adapter to the meter (Figure 2).
3. Connect the universal power adapter to an AC receptacle (Figure 3).
4. Turn the meter on.


Figure 2 Connector panel

- | |
|--------------------------------|
| 1 Magnetic stirrer connector |
| 2 Conductivity probe connector |

- | |
|----------------|
| 3 Power supply |
|----------------|


Figure 3 AC power connection

User interface and navigation

User interface

Keypad description


1 CALIBRATION key: start a calibration , view calibration data and change the calibration frequency	5 DOWN key: scroll to other options, change a value
2 TEMPERATURE and STIRRER key: change the temperature (°C, °F) or the stirring speed	6 ON/OFF: turn on or turn off the meter
3 UP key: scroll to other options, change a value	7 NaCl measurement key: start a salinity measurement
4 EC measurement key: start a conductivity measurement	

Display description

The meter display shows the concentration, units, temperature, calibration status, date and time.

Figure 4 Single screen display

1 Measurement unit and value (conductivity or salinity)	3 Sample temperature (°C or °F)
2 Measurement mode or time and date	4 Visual measurement timer

Navigation

Use the calibrate key to calibrate the probe. Use the parameter key to take a sample measurement. Use the TEMPERATURE and STIRRER key to change the temperature (°C, °F) or the stirring speed. Use the arrow keys to scroll to other options or to change a value. Be sure to look at the display during tasks as the screens change quickly. Refer to each task for specific instructions.

Start-up

Turn the meter on and off

NOTICE

Make sure that the probe is connected to the meter before the meter is turned on.

Push to turn on or turn off the meter. If the meter does not turn on, make sure that the AC power supply is properly connected to an electrical outlet.

Change the language

The display language is selected when the meter is powered on for the first time.

Use the ▲ or ▼ to select a language from a list. The meter returns automatically to the measurement screen after 3 seconds.

Note: To change the language from the main menu, push ▲.

Standard operation

Calibration




▲ WARNING

Chemical exposure hazard. Obey laboratory safety procedures and wear all of the personal protective equipment appropriate to the chemicals that are handled. Refer to the current material safety data sheets (MSDS) for safety protocols.

Calibration procedure

This procedure is for general use with liquid calibration solutions. Refer to the documents that are included with each probe for additional information.

Note: Solutions should be stirred during calibration. For more information about the stirring settings, refer to [Change the stirring settings](#) on page 9.

1. Pour the buffers or calibration solutions into the labeled calibration beakers.
2. From the main menu push  to select EC calibration.
3. Rinse the probe with deionized water and put the probe into the first calibration beaker. Be sure that there are no air bubbles in the probe measuring chamber.
4. Push  to measure the first calibration solution. The next calibration solution is shown.
5. Rinse the probe with deionized water and put the probe into the second calibration beaker. Be sure that there are no air bubbles in the probe measuring chamber.
6. Push  to measure the second calibration solution. The next calibration solution is shown.
7. Rinse the probe with deionized water and put the probe into the third calibration beaker. Be sure that there are no air bubbles in the probe measuring chamber.
8. Push  to measure the third calibration solution.

When the calibration is good, the display briefly shows 3rd Standard OK and then returns to the main menu.

View the calibration data

Data from the most recent calibration can be shown.

1. From the main menu push .
2. Use the ▼ to view the last calibration data. The meter returns automatically to the measurement screen after 3 seconds.

Set the calibration reminder

The calibration reminder can be set between 0 to 99 days (default 15 days). The display shows the remaining time to the new calibration.

Note: When 0 days is selected, the calibration reminder is turned off.

1. From the main menu push  and then push ▼ until the display shows Cal. frequency.
2. Use the ▲ or ▼ to change the value. The meter returns automatically to the measurement screen after 3 seconds.

Sample measurements

Each probe has specific preparation steps and procedures for taking sample measurements. For step-by-step instructions, refer to the documents that are included with the probe.

Note: Solutions should be stirred during measurement. For more information about the stirring settings, refer to [Change the stirring settings](#) on page 9.

Push NaCl or EC to take a sample measurement. During measurements, the parameter flashes and the timer shows the stabilization time. To change the measurement parameter (if applicable), push and hold NaCl or EC.

To measure the sample continuously, push NaCl or EC to start a measurement and push NaCl or EC again during the stabilization. The parameter flashes to indicate the continuous measurement mode.

Advanced operation

Change the date and time

The date and time can be changed from the Date and Time menu.

1. Use the ▼ to enter the Date and Time menu. The date and time format: dd-mm-yyyy 24h.
The attribute to be changed will be automatically highlighted and will advance to the next attribute if no ▲ or ▼ is pushed.
2. To change an attribute, push ▲ or ▼ when the highlighted attribute is selected.
The current date and time will be shown on the display.

Adjust the display contrast

1. Push the ▲ and ▼ at the same time to enter the Display contrast menu.
2. Use the ▲ or ▼ to adjust the contrast of the display. The meter returns automatically to the measurement screen after 3 seconds.

Change the stirring settings

The stirring speed can be changed during calibration and during a measurement.

1. Push $\frac{T}{S}$ during a calibration or during a measurement to enter the stirring menu.
2. Use the ▲ or ▼ to change the stirring speed in %.

Change the temperature units

The temperature units can be changed to Celsius or Fahrenheit.

1. From the main screen push $\frac{T}{S}$.
2. Push $\frac{T}{S}$ to select between Celsius or Fahrenheit.

Maintenance

▲ WARNING

Multiple hazards. Do not disassemble the instrument for maintenance or service. If the internal components must be cleaned or repaired, contact the manufacturer.

▲ CAUTION

Personal injury hazard. Only qualified personnel should conduct the tasks described in this section of the manual.

Clean the instrument

NOTICE

Never use cleaning agents such as turpentine, acetone or similar products to clean the instrument including the display and accessories.

Clean the exterior of the instrument with a moist cloth and a mild soap solution.

Clean the probe

Clean the probe as needed. Refer to [Troubleshooting](#) on page 11 for more information about cleaning. Refer to the probe documentation for information about the probe maintenance.



Use the cleaning agents listed in [Table 1](#) for contaminations on the conductivity probe.

Table 1 Cleaning agents for the conductivity probe

Contamination	Cleaning agent
Grease, oils, fats	Electrode cleaning solution
Limescale	0.1 N HCl solution

Replace the magnetic stirrer

If the magnetic stirrer does not start, follow the numbered steps to replace the magnetic stirrer.

Troubleshooting

Refer to the following table for common problem messages or symptoms, possible causes and corrective actions.

Table 2 Calibration warnings and errors

Error/Warning	Solution
UNSTABLE READING Time t > 100 s	Repeat with  Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol.
Check temperature. Check the cell.	Make sure the standard temperature is constant during the calibration. Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol.
STANDARD TEMPERATURE < 15 °C or > 35 °C	Repeat with  The standard temperature must be between 15 and 35 °C.
Difference C > 30 %	Repeat with 
C < 0.05 or > C > 50.00 °C	Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol; connect a different probe to verify if problem is with probe or meter. Examine the standard solution: Use new standard solution.
SAME BUFFERS	
Buffer or cell in poor conditions	
BUFFER NOT RECOGNIZED	
Difference C > 30 %	
SAME STANDARDS	

Table 3 Measurement warnings and errors

Error/Warning	Solution
23.2 °C 728 µS/cm @ 25 °C	Repeat with 
EC out of range	Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol; connect a different probe to verify if problem is with probe or meter.
Out of range °C	Examine the temperature sensor. Connect a different probe to verify if problem is with probe or meter.
TC = 0 Impossible to measure salinity	Modify TC
Time > 150 s	Examine the temperature. Examine the probe: Make sure that the probe is properly immersed in the sample; make sure that there are no air bubbles inside the probe measuring chamber; rinse the probe with ethanol; connect a different probe to verify if problem is with probe or meter.

Replacement parts and accessories

Note: Product and Article numbers may vary for some selling regions. Contact the appropriate distributor or refer to the company website for contact information.

Replacement parts

Description	Item no.
sensION+ PH3 Lab pH-meter with accessories, without probe	LPV2000.98.0002
sensION+ PH31 Lab pH-meter, GLP, with accessories, without probe	LPV2100.98.0002

Replacement parts and accessories (continued)

Description	Item no.
sensION+ MM340 Lab pH & Ion-meter, GLP, 2 channels, with accessories, without probe	LPV2200.98.0002
sensION+ EC7 Lab conductivity meter, with accessories, without probe	LPV3010.98.0002
sensION+ EC71 Lab conductivity meter, GLP, with accessories, without probe	LPV3110.98.0002
sensION+ MM374, 2 channel Lab meter, GLP, accessories, without probes	LPV4110.98.0002

Consumables

Description	Item no.
Conductivity standard solution 147 $\mu\text{S}/\text{cm}$, 125 mL	LZW9701.99
Conductivity standard solution 1413 $\mu\text{S}/\text{cm}$, 125 mL	LZW9711.99
Conductivity standard 12.88 mS/cm , 125 mL	LZW9721.99
Conductivity standard 147 $\mu\text{S}/\text{cm}$, 250 mL	LZW9700.99
Conductivity standard 1413 $\mu\text{S}/\text{cm}$, 250 mL	LZW9710.99
Conductivity standard 12.88 mS/cm , 250 mL	LZW9720.99
Enzyme solution	2964349
Pepsin Cleaning Solution	2964349
Electrode cleaning solution	2965249
0.1 N HCl solution	1481253

Accessories

Description	Item no.
Magnetic stirrer with sensor holder, for sensION+ MM benchtop	LZW9319.99
3x50 mL printed beaker for benchtop conductivity calibration	LZW9111.99
Three-sensor holder, for sensION+ benchtop instruments	LZW9321.99
Holder and clamp for three sensors	LZW9155.99
Pyrex glass chamber, continuous flow measurements	LZW9118.99
PP protector, electrode storage	LZW9161.99

Standard solutions

Conductivity standard solutions

Refer to [Table 4](#) for conductivity values of standard solutions at varying temperatures.

Table 4 Conductivity and temperature values

Temperature		Conductivity (EC)			
$^{\circ}\text{C}$	$^{\circ}\text{F}$	$\mu\text{S}/\text{cm}$	$\mu\text{S}/\text{cm}$	mS/cm	mS/cm
15.0	59	119	1147	10.48	92.5
16.0	60.8	122	1173	10.72	94.4
17.0	62.6	125	1199	10.95	96.3
18.0	64.4	127	1225	11.19	98.2
19.0	66.2	130	1251	11.43	100.1
20.0	68	133	1278	11.67	102.1
21.0	69.8	136	1305	11.91	104.0
22.0	71.6	139	1332	12.15	105.4
23.0	73.4	142	1359	12.39	107.9

Table 4 Conductivity and temperature values (continued)

Temperature		Conductivity (EC)			
°C	°F	µS/cm	µS/cm	mS/cm	mS/cm
24.0	75.2	145	1386	12.64	109.8
25.0	77	147	1413	12.88	111.8
26.0	78.8	150	1440	13.13	113.8
27.0	80.6	153	1467	13.37	115.7
28.0	82.4	156	1494	13.62	—
29.0	84.2	159	1522	13.87	—
30.0	86	162	1549	14.12	—
31.0	87.8	165	1581	14.37	—
32.0	89.6	168	1609	14.62	—
33.0	91.4	171	1638	14.88	—
34.0	93.2	174	1667	15.13	—
35.0	95	177	1696	15.39	—

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P.O. Box 389, Loveland, CO 80539-0389 U.S.A.
(800) 227-4224 (U.S.A. only)