## ESPAÑOL

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## ENGLISH

### Blood collection Eurotubo® tubes

Tubes are manufactured from virgin polypropylene. The tubes are inert to practically all chemicals associated with blood collection and are virtually unbreakable due to the excellent impact properties of polypropylene. The cap is manufactured from virgin polyethylene. The ergonomic design of the cap provides convenient one-handed cap removal. Lot number and Expiry date are coded on the tube label to provide 100% traceability.

#### Tubes for SERUM:

Serum tubes are provided with clot accelerator additives to reduce clotting time and are available with either gel or granular serum separators.

#### **Directions of use:**

- Remove the cap.
- Introduce the blood sample into the tube. Some assays, such as glucose, require specific blood volumes; fill these tubes up to the fill-line printed on the label.
- Recap the tube securely and blend the contents by gently turning the tube end-for-end four or six times.
- Set or rack the tube vertically and allow fifteen minutes to assure proper clotting.
- Centrifuge it at 1.500 xg for fifteen minutes.

#### Tubes for PLASMA:

Plasma tubes are provided with different additives specific to the analysis desired. The following additives are available: Lithium Heparine, Iodoacetate Heparin (ideal for biochemical tests), Citrate (for ESR and coagulation), EDTA Dipotassium and EDTA Tripotassium.

These tubes (except the citrate tube) include the additive in spray in order to help the mixture with the blood. Citrate tubes (ESR and coagulation) are manufactured with liquid anticoagulant according to the standard procedures.

#### **Directions of use:**

- Remove the cap.
- Introduce the blood into the tube and fill it up to the 'Fill-Line', printed on the label.
- Recap the tube securely and blend the contents by gently turning the tube end-for-end four or six times.
- Whole blood may be analysed immediately. For analysing plasma it is recommended to wait for some minutes or to spin the tube softly.

For more information refer to our general catalogue. Detailed information on each tube is available.

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