MSwab[®]

Instructions for Use

CE IVD

Copan MSwab® collection, transport and preservation system for molecular and culture applications

Instructions for use

INTENDED USE

The MSwab® system is used for the collection, transport and preservation of clinical specimens from the collection site to the testing laboratory. In the laboratory, specimens collected in MSwab® system can be analysed using standard clinical procedures for:

- bacterial culture of aerobic and facultative anaerobic gram-positive microorganisms;
- viral culture of HSV 1 and HSV 2 viruses;
- Nucleic acid detection of bacteria and viruses.

SUMMARY AND PRINCIPLES

One of the routine procedures in diagnosing bacterial infections involves the safe collection and transport of swab specimens. This can be accomplished by using the Copan MSwab® collection, transport and preservation system. The medium has been designed to maintain the viability of aerobic and facultative anaerobic gram-positive bacteria and HSV 1 and HSV 2 viruses during transport to the testing laboratory.

In addition, Copan MSwab® allows preserving nucleic acids of pathogens to be identified using molecular amplification techniques (NAAT).

Copan MSwab® Collection, Transport and Preservation System is supplied in three different formats:

- a) Collection kit format. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 1 ml or 3 ml of MSwab® transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwabs® that has a tip flocced with nylon fiber.
- b) Tube only format. A plastic screw-cap tube with conical shaped bottom filled with 1 ml, 2 ml or 3ml of MSwab® transport and preservation medium.
- c) Collection kit format with a cleaning swab. Each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 2 ml of MSwab® transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwabs® that has a tip flocced with nylon fiber and a sterile peel pouch containing one cleaning swab with large tip in fiber wound for the removal of excess vaginal mucus.

Once the specimen has been collected using the swab, it must immediately be inserted into the MSwab® test tube containing the transport and preservation medium. Specimens collected using MSwab®, to be analysed using bacterial or viral culture techniques, must be transported directly to the laboratory, preferably within 2 hours from collection^(1, 2, 7); in order to maintain optimal microorganisms viability. If immediate specimen delivery or analysis is delayed, specimens can be refrigerated at 4-8°C or stored at room temperature (20-25°C) and analysed within 48 hours from sampling. Bacterial viability studies for *Staphylococcus aureus*, ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (methicillin-resistant) ATCC® 43300 and ATCC® 700698 show viability of the tested microorganisms for up to 14 days in a refrigerated environment (4-8°C) or 72 hours at room temperature (20-25°C). Independent scientific studies on swab transport systems show that the viability of some bacteria is higher at refrigerated temperature than at room temperature⁽¹²⁻²¹⁾.

If viral samples must be frozen, they should be frozen at -70°C.

Specimens collected using MSwab® to detect bacterial or viral nucleic acids must be analysed within 14 days if stored at room temperature (20-25°C), within 21 days if stored at 4°C and within 6 months if stored at -20°C.

REAGENTS

MSwab® transport medium formulation

Distilled water
Organic solvent
Buffer
Bovine serum albumin

pH: 8.5 ± 0.20

REQUIRED MATERIALS BUT NOT PROVIDED

Materials suitable for aerobic and facultative anaerobic gram-positive bacteria isolation and culture.

These materials include culture plates or test tubes and incubation systems. Refer to the laboratory manuals for the recommended protocols relative to the techniques for aerobic and facultative anaerobic bacteria culture and identification from clinical specimen swabs^(2, 4). Materials suitable for virus isolation, differentiation and culture. These materials include cell lines for tissue culture, culture medium for tissues, incubation systems and reading instruments. Refer to the appropriate references for the recommended protocols for virus isolation and identification^(1, 7). Materials suitable for the extraction and amplification of nucleic acids for molecular biological testing. Heating and centrifuge block for quick extraction methods. Refer to laboratory reference manuals for identification and amplification of nucleic acids in clinical swab samples.

In the case of rapid extraction methods, it is suggested to transfer the appropriate aliquot of MSwab® (refer to the section **Processing MSwab® specimens for molecular testing in the laboratory, Method B**) to Copan test tubes with glass beads (code available separately 2E013S50).

PRODUCT STORAGE

The product is ready for use and does not require further preparation. It must be stored in its original container at 5-25°C until used. Do not overheat. Do not incubate or freeze prior to use. Incorrect storage will result in loss of effectiveness. Do not use after the expiry date, which is clearly printed on the outer box and on each individual collection unit and the specimen transport tube label.

SPECIMEN COLLECTION, TRANSPORT AND PRESERVATION

Specimens collected for microbiological investigations which comprise the isolation of bacteria or viruses should be collected and handled following published manuals and guidelines^(7, 6, 4).

To maintain optimum microorganisms viability and nucleic acids integrity, transport specimens collected using MSwab® directly to the laboratory, preferably within 2 hours of collection ^(1, 2, 7). If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 hours. Bacterial viability studies for *Staphylococcus aureus*, ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC® 43300 and ATCC® 700698 show viability of the tested microorganisms up to 14 days at refrigerated temperature (4 – 8°C) or 72 hours at room temperature (20 – 25°C). If viral samples must be frozen, they should be frozen at -70°C.

Specimens collected using MSwab® to be analyzed with NAAT should be processed within 14 days when stored at room temperature (20 – 25°C), within 21 days when stored at 4°C and within 6 months when stored at -20°C.

Performance testing with Copan MSwab® was conducted using laboratory strains suspensions spiked onto the FLOQSwabs® associated to the medium. Tested strains include: *B. pertussis* ATCC® 8467, *S. aureus* (methicillin resistant) ATCC® 43300, *E. coli* O157-H7 ATCC® 700728, *S. tiphimurium* ATCC® 14028, *C. difficile* ATCC® 9689, Influenza A virus ATCC® VR-822, Influenza B virus ATCC® VR 786, *Neisseria gonorrhoeae* ATCC® 43069 and *Chlamydia trachomatis* ATCC® VR880. Performance testing was not conducted using human specimens or human matrices.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations ^(34, 35, 36, 37). Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

MATERIAL SUPPLIED

Fifty (50) units of Copan MSwab® System are contained in a box and 6 boxes x 50 units are contained in a carton.

Copan MSwab® Collection, Transport and Preservation System is supplied in three different formats:

1. Collection kit format: each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 1 ml or 3 ml of MSwab® transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwabs® that has a tip flocced with nylon fiber (Fig.2a).
2. Tube only format: a plastic screw-cap tube with conical shaped bottom filled with 1 ml, 2 ml or 3 ml of MSwab® transport and preservation medium.
3. Collection kit format with a cleaning swab: each collection kit consists of a package containing a plastic screw-cap tube with conical shaped bottom filled with 2 ml of MSwab® transport and preservation medium and a small sterile peel pouch containing one specimen collection FLOQSwabs® that has a tip flocced with nylon fiber (Fig.2a) and a sterile peel pouch containing one cleaning swab with large tip in fiber wound for the removal of excess vaginal mucus (Fig.2b).

There are two types of specimen collection swabs: a standard-size swab with Nylon flocced tip intended for the collection of specimens from anatomical sites such as throat, vagina, wounds, rectum and faeces; a flexible minitip with Nylon flocced tip for nasopharyngeal collection. All flocced swabs supplied with MSwab® have a breakpoint on the swab shaft marked with a coloured line.


After collecting the specimen from the patient, the pre-moulded breakpoint facilitates breakage of the applicator into the test tube containing the MSwab® transport medium. The special shape of the capture caps and the test tubes allows capturing the swab shaft after it is broken off and the cap is closed (see Fig.1).

Screwing the cap onto the test tube, the end of the shaft is moved into the cavity of the cap. When the test tube is opened in the testing laboratory, the applicator remains attached to the cap and the operator can easily remove the swab from the test tube and perform microbiological analyses using the cap as a handle to hold the swab applicator.

The Capture Cap feature is only guaranteed when using the Copan standard-size flocced swab.

The Capture Cap feature is not applicable to the pernasal swab (REF. 6E013N and 6E092N01).

Fig. 1 - Capture of the broken swab shaft in the cap of the MSwab® test tube

LIMITATIONS

1. The conditions, time and volume of the specimens collected for culture are significant variables for obtaining reliable results. Follow the recommended guidelines for specimen collection ^(7, 8, 4).
2. MSwab® may not be used as enrichment, selective or differential medium.
3. MSwab® medium does not contain antibiotics. Patient specimens which may contain a high load of bacterial contaminants may require adding antibiotics to the cell preservation and culture medium.
4. The Copan MSwab® performance tests were conducted using laboratory strains applied on a swab following the test protocols described in Clinical and Laboratory Standards Institute M40-A2 Approved Standard ⁽⁹⁾. The performance tests were not conducted using human specimens.
5. The Copan MSwab® performance tests were conducted using Copan flocced swabs.
6. Take the necessary precautions against the biological risk and use approved aseptic techniques. The product may only be used by trained and qualified persons.

WARNINGS

1. When handling clinical specimens in the laboratory, wear gloves and any other protection devices necessary. When handling or analysing specimens from a patient, observe Biosafety Level 2 established by the CDC ^(51, 32, 33, 34).
2. For in vitro diagnostic use.
3. Do not re-sterilize unused swabs.

4. This product is for single use only; reuse may cause a risk of infection and/or inaccurate results.
5. Do not re-pack.
6. Not suitable for any other application than intended use.
7. Use of the product with a rapid diagnostic kit or with diagnostic instruments must be validated by the user a priori.
8. Do not use excessive force, pressure or bending when collecting swab samples from patients as this may result in breakage of the swab shaft.
9. Do not use the same test tube for more than one patient. This will result in incorrect diagnosis. Before transporting, make sure MSwab® screw cap is tightly closed.
10. Do not bend the collection swab before use.
11. Do not use the Mswab® medium to pre-moisten or pre-wet the swab applicator before collecting the specimen or to moisten the sampling sites.
12. Avoid contact of the MSwab® medium with skin and mucous membranes. If contact occurs, rinse with plenty of water.
13. Do not ingest the transport medium.
14. Carefully follow the instructions for use.
15. The product may only be handled by trained persons.
16. Observe approved biohazard precautions and aseptic techniques.
17. Repeated freezing and thawing of specimens may reduce the recovery of viable organisms ^(6,35).
18. Condition, timing and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.
19. It must always be assumed that all specimens contain infected microorganisms and the utmost caution is therefore recommended. After use, dispose of the test tubes and the swabs in accordance with the laboratory procedures for infected waste. Observe Biosafety Level 2 established by the CDC ^(31, 32, 33, 34). Discard unused reagents, waste and specimens according to local regulations.
20. Copan MSwab® must not be used if (1) there is evidence of product damage (e.g. broken swab tip or shaft), or contamination, (2) there is evidence of leakage, (3) the expiry date has passed, (4) the swab package is open, (5) there are other signs of deterioration.
21. Due to the geometry of the Flexible Minitip, the swab may coil when placed into the test tube, therefore, if should be necessary remove the swab from the tube, pay attention and observe appropriate biohazard precautions to protect the operator and the environment in case of splashes.
22. Check the version of the operating instructions. The correct version is the one supplied with the device or available in electronic format, and can be identified by the e-IFU indicator on the packaging label.

INSTRUCTIONS FOR USE

Copan MSwab® is ready for use and does not require further preparation. It is available in the different configurations shown in Table 1.

REF.	Copan MSwab® product description	Package
6E012N	Sample collection pack containing: - 12x80mm test tube with screw cap containing 1ml of MSwab® preservation and transport medium. - Standard-size FLOQSwabs® swab with Nylon flocked tip, sterile and individually wrapped.	50 units/box of 6x50 units/cardboard box
6E013N	Sample collection pack containing: - 12x80mm test tube with screw cap containing 1ml of MSwab® preservation and transport medium. - Flexible minitip with Nylon flocked tip, sterile and individually wrapped	50 units/box of 6x50 units/cardboard box
6E092N01	Single transport and preservation tube: - 16x100mm round-bottom test tube with screw cap containing 3 ml of MSwab® preservation and transport medium. - Flexible minitip swab with Nylon flocked tip, sterile and individually wrapped	50 units/box of 6x50 units/cardboard box
6E011N	Single transport and preservation tube: - 12x80mm tapered-bottom test tube with screw cap containing 1ml of MSwab® preservation and transport medium.	50 units/box of 6x50 units/cardboard box
6U019N	Single transport and preservation tube: - 12x80mm tapered-bottom test tube with screw cap containing 2ml of MSwab® preservation and transport medium.	50 units/box of 6x50 units/cardboard box
6E076N	Single transport and preservation tube: - 12x80mm tapered-bottom test tube with screw cap containing 3ml of MSwab® preservation and transport medium.	50 units/box of 6x50 units/cardboard box
6E028N.MER	Single use sample collection pack containing: - 12x80mm test tube with screw cap containing 2ml MSwab® storage and transport medium. - Standard-size FLOQSwabs® swab with Nylon flocked tip, sterile and individually wrapped. - Cleaning swab with Rayon fiber large tip, sterile and individually wrapped.	50 units/box of 6x50 units/cardboard box

Table 1

Specimen collection

Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious microorganisms. For specific guidance regarding specimen collection procedures, consult published reference manuals (7, 2).

For MSwab® codes 6E012N, 6E013N and 6E092N01:

1. Open the MSwab® kit package, remove the test tube containing the medium and the inner pouch containing the sterile swab.
2. Remove the sterile swab from its pouch (see Figure 2) and collect the clinical specimen. The operator must touch the swab applicator only above the colored breakpoint line, as illustrated in Figure 2, which is the opposite end to the nylon fiber tip. At all times when handling the swab applicator, the operator must not touch the area below the colored breakpoint line (the area from the line to the tip of the nylon flocked swab) as this will lead to contamination of the applicator shaft and the subsequent culture. To prevent the risk of contamination, ensure that the swab tip only comes into contact with the collection site.
3. Collect the sample from the patient.
4. Unscrew and remove the cap from MSwab® tube making sure not to spill the medium.
5. After collecting the specimen, insert the swab into the test tube until the breakpoint is at the level of the test tube opening. Bend the swab shaft at a 180 degree angle to break it at the red marked breakpoint. If necessary, gently turn the swab shaft until it is completely broken off and remove the upper part of the swab shaft (Figure 2).
6. Discard the broken off part of the swab shaft in an approved medical waste disposal container.
7. Replace cap on the tube and secure tightly.
8. Write the patient name and data on the test tube label or apply patient identification label.
9. Send the specimen to the laboratory.

For MSwab® code 6E028N.MER:

1. Open the kit package, remove the test tube containing the transport medium and the two inner pouches, one containing the sterile swab and one containing the "Cleaning Swab".
2. Use of the cleaning swab is recommended only in the case of an endocervical specimen collection procedure. The cleaning swab with large tip in fiber tip (see Figure 2.b) is only used to remove excess mucus present at the level of the uterine cervix opening and the surrounding mucous membranes. **After use, discard the swab in an approved sanitary waste container.** In case of procedures for collecting specimens other than endocervical specimens, the cleaning swab must not be used and must be discarded.
3. For clinical specimen collection, follow the instructions for codes 6E012N, 6E013N and 6E092N01 from step 2 onward.

For MSwab® codes 6E011N, 6U019N and 6E076N:

1. After collecting the specimen from the patient using a swab, unscrew and remove the cap from the MSwab® test tube being careful not to spill the transport medium.
2. Insert the swab into the test tube.
3. If the swab has a breakpoint, bend and break the swab at the breakpoint taking care to hold the test tube away from your face.
4. Replace cap on the tube and secure tightly.
5. Write the patient name and data on the test tube label or apply patient identification label.
6. Send the specimen to the laboratory.

Fig. 2 - Collection swab with breakpoint indication line and applicator handling area


Fig. 2.a COLLECTION SWAB



Fig 2.b CLEANING SWAB

Processing MSwab® specimens in the laboratory – Bacteriology

MSwab® samples should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines⁽¹⁻⁶⁾.

Culture investigations of swab specimens for the presence of bacteria routinely involve the use of solid agar culture medium in Petri dish plates. The procedure for inoculation of MSwab® samples onto solid agar in Petri dishes is as follows.

Note: When handling clinical specimens, wear latex gloves and any other protection devices necessary. Observe the other Biosafety Level 2 recommendations issued by the CDC^(31, 32, 33, 34).

Collection Kit Format (with or without cleaning swab):

Vortex the MSwab® test tube containing the specimen collected with the swab for 5 seconds and evenly disperse and suspend the patient specimen in the transport medium.

1. Unscrew the MSwab® cap and remove the swab applicator.
2. Roll the tip of the MSwab® applicator onto the surface of one quadrant of the culture media plate to provide the primary inoculum.
3. If it is necessary to culture the swab specimen onto a second culture media plate, return the MSwab® applicator to the transport medium tube for two seconds to absorb and recharge the applicator tip with transport medium/patient sample suspension then repeat Step No. 3.
4. If it is necessary to inoculate additional culture media plates, return the MSwab® applicator to the transport medium tube and recharge the swab applicator tip with the transport medium/patient sample suspension before inoculating each additional plate.

The procedure described above utilizes the MSwab® applicator like an inoculation wand to transfer the suspension of patient sample in transport medium onto the surface of a culture plate creating the primary inoculum (see Fig 3).

Alternatively, the operator can vortex mix the MSwab® tube with the swab inside for 5 seconds and then transfer 100µl volumes of the suspension onto each culture plate using a volumetric pipettor and sterile pipet tips. Standard laboratory techniques should then be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig 4).


Tube Only Format:

Vortex the MSwab® test tube containing the specimen for 5 seconds and evenly disperse and suspend the patient specimen in the transport medium. Transfer a minimum of 100µl volumes of the suspension onto each culture plate using a volumetric pipettor and sterile pipet tips. Standard laboratory techniques should then be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig 4).

Fig. 3 - Inoculation procedures of MSwab® specimens on solid agar in Petri dishes


1. Use the swab to inoculate the specimen.

2. Use a pipettor and sterile pipet tips to inoculate a minimum of 100µl of specimen.

This inoculation procedure should be considered preferable if the clinical specimen is also analysed using NAAT techniques.

Fig. 4 - Procedure for streaking MSwab® specimens on agar Petri dishes for primary isolation⁽³³⁾

Seed a primary inoculum of MSwab® specimen onto the surface of an appropriate agar culture plate in the first quadrant.

Use a sterile bacteriology loop to streak the primary inoculum across the surface of the second, third and fourth quadrants of the agar culture plate.

Preparation of Gram Stain Smears of MSwab® Specimens

Laboratory analysis of clinical swab samples collected from certain sites on the patient can routinely include microscopic examination of stained preparations ("direct smears") using the Gram stain procedure. This can provide valuable information to physicians who are managing patients with infectious diseases⁽²²⁾. There are many instances in which a Gram stain can assist in making a diagnosis^(23, 27).

The Gram stain can also help to judge specimen quality and contribute to the selection of culture media especially with mixed flora.

Microscope slides of patient specimens transported in Copan MSwab[®] transport system can be prepared for Gram stain analysis, as describe below, by sampling an aliquot of vortexed suspension of the swab^(3,4). Sample transported in MSwab[®] elution medium represent a homogeneous suspension in liquid phase. It can be uniformly smeared allowing clear and easy reading.

Note: When handling clinical specimens, wear latex gloves and any other protection devices necessary. Observe the other Biosafety Level 2 recommendations issued by the CDC^(31, 32, 33, 34).

1. Take a clean glass microscope slide, place it on a flat surface and inscribe an area using a diamond-tipped or similar glass marker to identify the location of the specimen inoculum. Note: a slide with a pre-marked 20 mm well can be used.
2. Vortex mix the MSwab[®] tube containing the swab sample for 5 seconds to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the medium.
3. Unscrew the MSwab[®] cap and using a sterile pipet, transfer 1 – 2 drops of sample suspension to the inscribed area on the glass slide. Note: about 30µl would be a suitable amount of liquid for a pre-marked 20 mm diameter well slide.

In case of bloody or thicker specimens particular care should be taken to thinly spread the sample on the slide. Bacteria are difficult to detect if the sample shows many red cells and debris.

4. Allow the specimen on the slide to air dry at room temperature or place the slide in an electric slide warmer or incubator set at a temperature not exceeding 42°C.
5. Fix smears using methanol. Methanol fixation is recommended as it prevents lysis of Red Blood Cells, avoids damage to all host cells and results in a cleaner background^(3, 4, 22).
6. Follow published laboratory reference manuals and guidelines for performing the Gram stain. If commercial Gram stain reagents are used, it is important to comply with instructions in the manufacturer's product insert for performance test procedure.

For further information or guidance on the preparation of specimen slides for microscopic analysis, for information on Gram staining procedures and the interpretation and reporting of microscopic analysis, consult published laboratory reference manuals^(1 - 5, 22 - 27).

Processing MSwab[®] specimens in the laboratory – Virology

The survival of HSV 1 and HSV 2 depends on many factors including the type and concentration of the microorganism, duration of transport and storage temperature. To maintain optimum viability, specimens should be transported directly to the laboratory preferably within 2 hours of collection^(1, 2, 7, 29). If immediate delivery or processing is delayed, then specimens collected using Copan MSwab[®] Collection, Transport and Preservation System should be refrigerated at 4-8°C or stored at room temperature (20-25°C) and processed within 48 hours. If samples must be frozen, they should be frozen placed at -70°C.

In simulated transport and storage studies, the Copan MSwab[®] System was shown to be capable of maintaining the viability of HSV 1 and HSV 2 at refrigerated (4-8°C) and room temperature (20-25°C) conditions for up to 48 hours. Based on performance studies conducted by Copan and independent scientific publications, viability of certain microorganisms is superior at refrigerated temperatures compared with room temperature^(12 - 21, 29).

MSwab[®] specimens should be processed for virology culture using recommended cells line and laboratory techniques that will depend on the specimen type and the microorganism under investigation. For recommended shell vials and techniques for the isolation and identification of HSV 1 and HSV 2 from clinical swab specimens, refer to published virology manuals and guidelines^(1 - 6, 29, 30).

Culture investigations of swab specimens for the presence of HSV 1 and HSV 2 routinely involve the use of cells culture in shell vials. The procedure for inoculation of MSwab[®] specimens onto shell vials is described below.

NOTE: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens. Observe other BSL 2 recommendations.

Collection Kit Format:

1. Shake using a vortex mixer for 5 seconds the MSwab[®] tube containing the swab specimen to release the specimen from the swab tip and evenly disperse and suspend the patient specimen in the liquid transport medium.
2. Unscrew the MSwab[®] cap and remove the swab applicator.
3. Transfer 200 µl volumes of the suspension into the shell vial and proceed as per laboratory internal procedure.
NOTE: Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.
4. Proceed with the appropriate techniques for virus detection.

Tube Only Format:

1. Shake using a vortex mixer for 5 seconds the MSwab[®] tube.
2. Unscrew the MSwab[®] cap and transfer 200 µl volumes of the suspension into the shell vial and proceed as per laboratory internal procedure.
NOTE: Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.
3. Proceed with the appropriate techniques for virus detection.

Processing MSwab[®] specimens for molecular analysis in the laboratory

Specimens received in the laboratory for nucleic acid detection should be processed when received in the laboratory. In case of delay, please refer to the appropriate specimen storage conditions.

NOTE: Use latex gloves and other general means of protection to handle clinical specimens. Follow Biosafety Level (BSL) 2 established by the CDC^(31, 32, 33, 34).

When working with molecular methods care should be taken to prevent carry over contamination. Spatial separation of working areas and unidirectional workflow are essential to prevent amplicon carry-over⁽³⁵⁾.

MSwab[®] medium formulation has been developed to be compatible with PCR Master Mixes; this makes it suitable for direct NAAT assays without the need of a standard extraction and purification step.

Two different methods (A and B) can be applied to MSwab® sample processing:

A) Standard extraction method

1. Mix the MSwab® test tube in a vortexer for 10 seconds;
2. Unscrew the cap (NOTE: for tube only format go to point 4) and, holding it between thumb and index finger, spin it to extract most of the fluid from the tip.
3. Discard the swab;
4. Transfer the appropriate amount of specimen (e.g. 200 µl) to an extraction tube following normal laboratory procedures. For code 6E013N and 6E092N01 where the capture cap feature is not available, use tweezers to extract the applicator from the tube. Use caution and observe adequate biohazard precaution to protect the operator and the environment in case of splash.
5. Continue in accordance with the extraction and amplification procedures of the kits used.

The MSwab® system has been validated using the following extraction methods: silica-gel membrane (Qiagen, Macherey & Nagel) and magnetic beads (easyMAG, Qiasymphony, Magnapure). Other extraction methods are applicable after validation by the end user.

B) Rapid extraction method

A thermal shock step is suggested to lyse samples for viral detection.

1. Vortex the original test tube of the MSwab® specimen for 10 seconds.
2. Using a micropipette, transfer 200µl of the MSwab® specimen to a sterile micro test tube containing glass beads (2E013S50: Copan test tube with glass beads, available separately) and vortex for 10 seconds. Store the original MSwab® specimen for culture or additional testing.
3. Heat the micro test tube in a thermostat set to 98-100°C for at least:
 - 3 min FOR VIRUSES
 - 10 min FOR BACTERIA GRAM POSITIVE
 - 5 min FOR BACTERIA GRAM NEGATIVE
4. Cool the test tube at room temperature for 5 minutes
5. Centrifuge the micro test tube at 10000xg for 2 minutes to settle cell debris.
6. Transfer an aliquot of the processed MSwab® specimen to the master mix and proceed with the amplification step in accordance with the amplification kit procedure.

Rapid extraction method as described above has been tested with R-Biopharm RIDA®GENE Real Time PCR kits, with Seegene Anyplex FluA/B Typing Real Time detection kit and with Genesig Real Time PCR kits by Primer Design in comparison with standard extraction method. Other PCR amplification kit may be also applicable prior validation. Tested strains include: B. pertussis (ATCC® 8467), S. aureus (methicillin resistant) (ATCC® 43300), E. coli O157-H7 (ATCC® 700728), S. typhimurium (ATCC® 14028), C. difficile (ATCC® 9689), Influenza A virus (ATCC® VR-822), Influenza B virus (ATCC® VR 786), Neisseria gonorrhoeae (ATCC® 43069) and Chlamydia trachomatis (ATCC® VR880). Tests were not conducted using human matrices. Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

QUALITY CONTROL

MSwab® applicators are tested to ensure they are non-toxic to bacteria. MSwab® medium and applicators are tested to ensure they are non-toxic to cell lines used for HSV 1 and HSV 2 cultivation. MSwab® transport medium is tested for pH stability and for absence of inhibitors with direct nucleic acid amplification⁽⁹⁾. MSwab® is quality control tested before release for ability to maintain viable Gram positive aerobic and facultative anaerobic bacteria and HSV viruses at room temperature (20 – 25°C) for specified time points.

Procedures for quality control of microbiology transport devices should be conducted using testing methods described in Clinical and Laboratory Standards Institute M40-A2 and other publications⁽⁹⁾. If aberrant quality control results are noted, patient results should not be reported.

RESULTS

The results obtained will depend largely on proper and adequate specimen collection as well as prompt transport and proper processing in the laboratory.

PERFORMANCE CHARACTERISTICS

Bacterial recovery

The test procedures employed for determining bacterial viability performance were based upon the quality control methods described in Clinical and Laboratory Standards Institute M40-A2⁽⁹⁾.

MSwab® system has an intended use limited to Gram positive aerobic and facultative anaerobic bacteria and HSV 1 and HSV 2, therefore its field application is more restricted than some other devices. For this reason the bacterial recovery studies were conducted under the simulated transport and storage conditions as described and defined in CLSI M40-A2, Quality Control of Microbiological Transport Systems: Approved Standard and included only the Gram positive aerobic and facultative anaerobic strains from the Group 1 of paragraph 7.11.1 of the CLSI M40-A2 document, in particular:

Streptococcus pyogenes	ATCC® 19615
Streptococcus pneumoniae	ATCC® 6305
Pseudomonas aeruginosa	ATCC® BAA-427

In addition, Copan included testing of additional Gram positive aerobic and facultative anaerobic microorganisms, clinically relevant, not required by CLSI M40-A2. The specific bacterial strains used in these studies are here listed:

Enterococcus faecalis	ATCC® 29212
Staphylococcus epidermidis	ATCC® 12228
Staphylococcus aureus	ATCC® 29213
Streptococcus agalactiae (Group B Strep)	ATCC® 13813
Listeria monocytogenes	ATCC® 19114
Bacillus cereus	ATCC® 10876
Staphylococcus aureus (Methicillin resistant)	ATCC® 43300
Staphylococcus aureus	ATCC® 6538
Staphylococcus aureus (Methicillin resistant)	ATCC® 700698

All cultures of bacteria were ATCC® (American Type Culture Collection) and were obtained commercially.

The selection of these organisms also reflects those Gram positive aerobic and facultative anaerobic bacteria that would normally be encountered in specimens collected and analyzed in a typical clinical microbiology laboratory.

Bacterial viability studies were performed on the Copan MSwab® at two different ranges of temperature, 4 – 8°C and 20 – 25°C, corresponding to refrigerator and room temperature, respectively. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs, 24 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate or Swab Elution Method.

Additional bacterial viability studies for *Staphylococcus aureus*, ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC® 43300 and ATCC® 700698 were performed on the Copan MSwab® at two different ranges of temperature, 4 – 8°C and 20 – 25°C, corresponding to refrigerator and room temperature, respectively.

Swab accompanying the transport system was inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and:

For studies performed at 4 – 8°C inoculated MSwab® tubes were held for 0 hrs, 10 days and 14 days. At the appropriate time intervals, each MSwab® was processed according to the Roll-Plate Method.

For studies performed at 20 – 25°C inoculated MSwab® tubes were held for 0 hrs and 72 hrs. At the appropriate time intervals, each MSwab® was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed on the Copan MSwab® at 4 – 8°C, corresponding to refrigerator temperature. Swabs accompanying each transport system were inoculated in triplicate with 100µl of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 hrs and 48 hrs. At the appropriate time intervals, each swab was processed according to the Roll-Plate Method.

Bacterial overgrowth studies were performed using *Pseudomonas aeruginosa*.

Viral infectivity studies were performed using HSV 1 and HSV 2. Swabs accompanying each transport system were directly inoculated in triplicate with 100µl of viral suspension. Swabs were then placed in their respective transport medium tubes and were held for 0, 24 and 48 hours at both 4°C and room temperature (20-25°C). At the appropriate time interval, each swab was vortexed, removed from its transport medium tube and then 200µl aliquots of this suspension was inoculated into shell vials. All cultures were processed by standard laboratory culture technique and examined after a specified incubation time. Viral infectivity was determined by fluorescing foci counts.

Viruses evaluated were:

Herpes Simplex Virus Type 1 (HSV 1)	ATCC® VR-539
Herpes Simplex Virus Type 2 (HSV 2)	ATCC® VR-734.

In accordance with Clinical and Laboratory Standards Institute M40-A2, viability performance is measured for each test microorganism at the 48 hrs time point and compared with the acceptance criteria.

In both the Roll-Plate and Swab Elution viability performance studies, Copan MSwab® System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) and room temperature (20 – 25°C). Acceptable recovery for the Roll-Plate Method is defined as ≥5 CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU. Acceptable recovery for the Swab Elution Method is defined as no more than a 3 log₁₀ (1 x 10³ +/- 10%) decline in CFU between the zero-time CFU count and the CFU of the swabs after the specified holding time.

Additional time-points were tested for *Staphylococcus aureus*, ATCC® 29213 and ATCC® 6538, and *Staphylococcus aureus* (Methicillin resistant) ATCC® 43300 and ATCC® 700698.

In the Roll-Plate viability performance studies, Copan MSwab® System was able to maintain acceptable recovery of all microorganisms evaluated at both refrigerator (4 – 8°C) for 14 days and room temperature (20 – 25°C) for 72 hrs. Acceptable recovery for the Roll-Plate Method is defined as ≥5 CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU.

Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8°C). For the Swab Elution Method, an overgrowth assessment is made on all bacteria species tested at the 48 hrs holding time point. Overgrowth assessment using the Swab Elution Method is defined as greater than 1 log₁₀ increase in CFU between the zero-time CFU count and the holding time point. For the Roll-Plate Method, an overgrowth assessment is made with a separate analysis in which swabs are dosed with 100µl containing 10² CFU of *Pseudomonas aeruginosa* culture. Overgrowth under these conditions is defined as greater than 1 log₁₀ increase in CFU between zero-time CFU and the 48 hrs holding time point.

Copan MSwab® System demonstrated no overgrowth based on the acceptance criteria described in Clinical and Laboratory Standards Institute M40-A2.

Copan MSwab® System was able to maintain the infectivity of the following viruses for at least 48 hours at both room temperature (20-25°C) and in the refrigerator (2-8°C) under the test conditions described above: Herpes Simplex Virus Type 1, Herpes Simplex Virus Type 2.

Nucleic acid preservation

Specimens collected using MSwab® to be analysed using nucleic acid amplification techniques must be processed within 14 days of sampling when stored at room temperature (20-25°C), within 21 days when stored at 4°C and within 6 months if stored at -20°C.

Performance tests with Copan MSwab® were conducted using suspensions of laboratory strains placed on the FLOQSwabs® swab associated with the medium. Tested strains include: *B. pertussis* ATCC® 8467, *S. aureus* (methicillin-resistant) ATCC® 43300, *E. coli* O157-H7 ATCC® 700728, *S. tiphimurium* ATCC® 14028, *C. difficile* ATCC® 9689, Influenza A virus ATCC® VR-822, Influenza B virus ATCC® VR 786, *N. gonorrhoeae* ATCC® 43069 and *C. trachomatis* ATCC® VR880. Performance tests were not conducted using human clinical specimens and/or matrices of human origin.

The results obtained mainly depend on appropriate and adequate specimen collection, as well as prompt transport and processing in the laboratory.

TABLE OF SYMBOLS

See the table of symbols at the end of the instructions for use.

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BIBLIOGRAPHY

1. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. *Manual of Clinical Microbiology*. 7th edition. Washington, DC: ASM; 1999.
2. Isenberg HD, Schoenkencht FD, Von Graevenitz A, Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, S. J. Rubin. American Society for Microbiology, Washington, DC, 1979.
3. Isenberg HD. (Editor in Chief). *Clinical Microbiology Procedures Handbook*. American Society for Microbiology, Washington, DC, 1992.
4. Isenberg HD. (Editor in Chief). *Essential Procedures for Clinical Microbiology*. American Society for Microbiology, Washington DC, 1998.
5. Summanen P, Baron EJ, Citron D, Strong C, Wexler HM, Finegold SM. (1993). *Wadsworth Anaerobic Bacteriology Manual*, 5th edn. Star Publishing Company, Belmont, CA.
6. ASM Cumitech publications number 10A, 13A, 17A. American Society for Microbiology, Washington, DC.
7. Miller JM. *A Guide to Specimen Management in Clinical Microbiology*. Second Edition. American Society for Microbiology. Washington, DC. 1999.
8. Miller JM, Holmes HT. Specimen collection, transport, and storage. In: *Manual of Clinical Microbiology*. 11th ed. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover FC, eds. Washington, DC: ASM; 1995:19-20.
9. Clinical and Laboratory Standard Institute (CLSI), 2014. *M40-A2 Quality Control of Microbiology Transport Systems; Approved Standard-Second Edition*
10. Sng E-H, Rajan VS, Teo K-L, Goh A-J. The recovery of *Neisseria gonorrhoeae* from clinical specimens: effects of different temperatures, transport times, and media. *Sex Trans Dis*. 1982; 9:74-78.
11. Sun Y, Taylor T, Williams L, Sautter RL. Comparison of bacterial viability using both the EZ brand collection and transport system with the Difco swab transport pack. Presented at: 96th ASM General Meeting. 1996; Washington DC. Abstract C35.
12. Arbique JC, Forward KR, LeBlanc J. Evaluation of four commercial transport media for the survival of *Neisseria gonorrhoeae*. *Diagnostic Microbiology and Infectious Disease*. 2000; 36:163-168.
13. Perry JL. Effects of temperature on fastidious organism viability during swab transport. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-55.
14. Wilson DA, Tuohy MS, Procop GW, Hall GS. Effects of storage on the recovery of bacteria from three swab transport systems: BD CultureSwab, BD Culturette and Starplex StarSwab II. 101st General Meeting of the American Society for Microbiology. 2001; Orlando, FL. Abstract C-61.
15. Arbique J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies described in NCCLS document M40-P Quality Control of Microbiology Transport Devices. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-40.
16. Mitchell E, Berman M, Ginocchio CC. Evaluation of two new Liquid Stuart transport systems: Platinum StarSwab II (Starplex Scientific) and BBL CultureSwab (Becton Dickinson). 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-74.
17. Perry JL, Matthews JS. Compliance of two popular swab transport systems with performance standards detailed by the new NCCLS Proposed Standard, M40-P. 103rd General Meeting of the American Society for Microbiology. 2003; Washington, DC. Abstract C-42.
18. Robinson A, Gruver ML. Comparison of bacterial survival in two transport systems stored at room temperature and refrigerator temperatures. 102nd General Meeting of the American Society for Microbiology. 2002; Salt Lake City, UT. Abstract C-69.
19. Human RP, Jones GA. Evaluation of 4 transport systems against a published standard. 104th General Meeting of the American Society for Microbiology. 2004; New Orleans, LA. Abstract C-161.
20. Human RP, Jones GA. Evaluation of swab transport systems against a published standard. *J Clin Pathol* 2004; 57:762-763.
21. Arbique J, Campbell S, MacFarlane M, Davidson RJ. Comparison of methodologies for anaerobic organisms described in NCCLS document M40-P, Quality Control of Microbiology Transport Devices. 13th European Congress of Clinical Microbiology and Infectious Disease (ECCMID). 2003; Glasgow, UK. Abstract P-652.
22. Marler LM, Siders JA, Allen SD. *Direct Smear Atlas, A Monograph of Gram-Stained Preparations of Clinical Specimens*. Lippincott Williams and Wilkins, 2001.
23. Rotimi VO, Yakubu Z, Abudu OO, Banjo TO. Direct Gram's stain of vaginal discharge as a means of diagnosing bacterial vaginosis. *Journal of Medical Microbiology*, 1991 Vol 35, Issue 2 103-106.
24. Spiegel CA, Amsel R, Holmes KK. Diagnosis of bacterial vaginosis by direct gram stain of vaginal fluid. *J. Clin Microbiol*. 1983 Jul; 18 (1):170-177.
25. Benavides MI, Moncada X, Rodriguez B, Castillo C. Gonococcal urethritis in men: clinical experience in 1978-1988. *Rev Med Chil*. 1992 Oct; 120(10):1140-3.
26. Mayaud P, Msuya W, Todd J, Kaatano G, West B, Begkoyian G, Grosskurth H, Mabey D. Rapid assessment in Rwandan refugee camps in Tanzania. *Genitourin Med*, 1997 Feb; 73 (1):33-8
27. Lynne S. *Garcia Clinical Microbiology Procedures Handbook, Third Edition 2010* ASM Press 3.2.1
28. AE Greenberg, LS Clesceri, and AD Eaton. 9215 heterotrophic plate count. In: *Standard Methods for the Examination of Water and Waste Water*. 18th ed. Washington, DC APHA; 1992: 9-33-9-34
29. Henry D. *Isenberg Essential procedures for Clinical Microbiology 1998* ASM PRESS 477-478
30. Steven Specter, Richard L. Hodinka, Stephen A. Young *Clinical Virology Manual Third edition*, 384-409
31. Fleming D. *Biological Safety: Principles and Practices*. January 2000. ASM, Washington DC.
32. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA. <http://www.cdc.gov/od/ohs/symp5/jrtext.htm>.
33. Richardson JH. *Biosafety in Microbiological and Biomedical Laboratories*. December 1994. Diane Publishing Company.
34. Hansen DJ. *Healthcare, Laboratories and Biosafety*. Vol 2., 1992. CRC Press.
35. *A Guide to Utilization of the Microbiology Laboratory for Diagnosis of Infectious Diseases: 2013 Recommendations by the Infectious Diseases Society of America (IDSA) and the American Society for Microbiology (ASM)*.

INDEX OF SYMBOLS

Symbol / Simbolo / Símbolo / Symbol / Symbole / Símbolos / Símbol / Символ	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Značenje / Значення
	Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Proizvođač / Виробник
	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / Dijagnostički uređaj u erpuveti / Медичний виріб для діагностики in vitro
	CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE/ Marcação CE / Oznaka CE / Маркування CE
	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Identifikasjon for godkjenningsorganet / Ідентифікаційний номер уповноваженого органу
	Sterilized using ethylene oxide / Sterilizzato usando ossido di etilene / Esterilizado por óxido de etileno / Sterilisiert mit Äthylenoxid / Stérilisé à l'aide d'oxyde d'éthylène / Esterilizado utilizando óxido de etileno / Sterilizovano korišćenjem etilen oksida / Стерилізований завдяки етиленоксиду
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Nemojte ponovo koristiti / Не використовувати повторно
	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Número do catálogo / Kataloški broj / Номер за каталогом
	Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Limites de temperature / Limites de temperatura / Opseg temperature / Температурні обмеження
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Upotrebiti do / Використати до
	Peel / Strappare per aprire / Desprender / Abziehen / Décoller / Destacar / Ogulite / Відкривати тут
	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Serijski (lot) broj / Код серії (партії)
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contém o suficiente para <n> testes / Sadržaj dovoljan za testova / Придатий для проведення <n> кількості випробувань
	Do not use if package is damaged / Non utilizzare in caso di confezionamento danneggiato / No utilizar en caso de paquete dañado / Bei Beschädigung der Verpackung nicht verwenden / Ne pas utiliser si l'emballage est abîmé / Não utilizar se a embalagem estiver danificada / Nemojte koristiti ako je ambalaža oštećena / Не використовувати, якщо упаковка пошкоджена
	This only applies to US: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Vale solo per gli Stati Uniti: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Sólo se aplica a los EE.UU.: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Gilt nur für die USA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / S'applique uniquement aux États-Unis: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Isto aplica-se apenas aos EUA: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / Ovo se odnosi samo na SAD: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner." / "Це стосується винятково США: "Caution: Federal Law restricts this device to sale by or on the order of a licensed practitioner."
	Date of Manufacture / Data di fabbricazione/ Fecha de fabricación / Herstellungsdatum/ Date de fabrication / Data de fabrico / Datum proizvodnje / Дата виробництва

Consult the operating instructions supplied with the device or available in electronic format, and which can be identified by the e-IFU indicator on the packaging label / Consultare le istruzioni per l'uso fornite con il dispositivo oppure disponibili in formato elettronico ed identificate dall'e-IFU indicator sull'etichetta imballo / Consultar las instrucciones de uso suministradas con el dispositivo o disponibles en formato electrónico e identificadas por el indicador e-IFU de la etiqueta del embalaje / Konsultieren Sie die Gebrauchsanweisung, die mit dem Gerät geliefert wird oder in elektronischem Format vorliegt und durch den e-IFU-Indikator auf dem Verpackungsetikett gekennzeichnet ist / Voir le mode d'emploi fourni avec l'appareil ou disponible au format électronique et identifiable grâce à l'indicateur e-IFU sur l'étiquette de l'emballage / Ver as instruções de utilização fornecidas com o dispositivo ou disponíveis em formato eletrônico e identificadas pelo indicador e-IFU na etiqueta da embalagem / Pogledajte uputstva za upotrebu koja ste dobili uz uređaj ili su dostupna u elektronskom formatu a koja se mogu prepoznati pomoću indikatora e-IFU na nalepnici na pakovanju / Див. друковані інструкції з використання або інструкції з використання в електронному форматі, якщо наявна позначка "eIFU Indicator"

Copan

Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia, Italy
Tel +39 030 2687211
Fax +39 030 2687250

Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:

Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com