Transsystem™

HPC072 Rev. 00 Date 2016.06

Transsystem™

Product Insert & How to Use Guide

**COPAN VIRAL TRANSYSTEM
COLLECTION AND TRANSPORT SYSTEM
Information Leaflet & User Guide**

See the glossary of symbols at the end of the information leaflet

INTENDED USE

Copan Viral Transystem is a sterile ready-to-use system for the collection, transport and preservation of clinical specimens for viral culture.

SUMMARY AND PRINCIPLES

One of the routine procedures in the diagnosis of infections caused by viruses, involves the collection and safe transportation of a clinical swab specimen from the patient to the laboratory. This can be accomplished using the Copan Viral Transystem. Each Copan Viral Transystem unit comprises a sterile peel pouch containing a rayon-tipped or polyester-tipped swab used to collect the sample and a tube containing a polyurethane foam sponge soaked with liquid transport medium. The foam sponge functions as a simple holding reservoir for the transport medium. The function of the sponge soaked medium is to keep the swab moist during transport. After sampling, the swab applicator is placed inside the tube and the tip of the swab makes contact with the sponge reservoir.

Copan Viral Transystem medium consists of a balanced salt solution for maintaining osmotic pressure within physiological limits and phosphate buffers to stabilize the pH of the medium. Viral Transystem is available with different applicator size and shaft, which facilitate the collection of specimens from various sites of the patient, as described below in the Directions for Use section. For specific recommendations about the collection of specimens for viruses and primary isolation techniques, consult the following ASM publications:
Lynne S. Garcia, Clinical Microbiology Procedures Handbook, Third Edition 2010 ASM Press 3.2.1

Once a swab sample is collected, it must be placed immediately into the transport tube where it comes into contact with the transport medium. Swab specimens for virus isolation and/or detection must be submitted to the laboratory as quickly as possible after collection. If immediate delivery is not possible, refrigerate specimens (2 - 8°C), or place them on wet ice or in a cold pack.

PRODUCT DESCRIPTION

250 units of Sterile Viral Transystem are contained in each box. Each individual package contains an applicator and a plastic tube containing the medium-soaked sponge.

Cat No.	Applicator Swab Type	Suggested Sampling Sites*
147C	Plastic/polyester Tip	Throat, cervical, vulvar, vaginal, rectal and nasal Dermal lesions and lesions of genital or mouth mucosa
148C	Aluminum Wire/polyester Tip	Eye/conjunctival, urethral
149C	Flexible Twisted Wire/ viscose Tip	Nasopharyngeal

* This is just a suggested table. Performance testing with Viral Transystem was conducted using laboratory strains spiked onto the swab. Performance testing was not conducted using human specimens. Please refer to your internal procedures to choose the most appropriate device for the specific sampling site.

REAGENTS

Copan Viral Transport medium incorporates the following components:

- Buffering agents
- HBSS
- Osmotized water

REQUIRED MATERIALS THAT ARE NOT INCLUDED

Appropriate materials for isolating, differentiating and culturing viruses. These materials include tissue culture cell lines, tissue culture medium, antibiotics, incubation systems and reading equipment.

STORAGE OF THE PRODUCT

Store Viral Transystem at 5-25°C. Do not freeze or overheat. Do not use after expiration date, which is clearly printed on the outer box and on each individual sterile pouch unit and on the specimen transport tube label.

LIMITATIONS

Viral Transystem is intended for the collection and transport of virology swab samples only. Samples for bacterial or Chlamydia investigation should be transported using alternative specific transport systems.

WARNINGS AND PRECAUTIONS

1. For In Vitro Diagnostic Use.
2. Directions for use must be followed carefully.
3. The swab applicator is qualified for short transient contact with the patient in order to collect a specimen. This short contact is made with the patient's external surfaces or internal surfaces via normal body orifices such as the nose, throat or vagina or via surgical wounds.
4. When collecting swab samples from patients, care should be taken not to use excessive force or pressure which may result in breakage of the swab shaft.
5. The fiber attachment to the applicator stick is qualified to withstand short transient contact with the patient in order to collect the sample; prolonged contact must be avoided as this may result in the detachment of the fiber.

6. When the swab sample is cultured in the laboratory, if the procedure necessitates the applicator(s) to be placed in a tube of culture broth, special attention should be taken in detaching the applicator stick from the cap to eliminate any risk of splashes or aerosols. If it is necessary to cut the applicator stick, sterile scissors should be used to facilitate a safe and clean break.
7. It must be assumed that all the specimens contain infectious micro-organisms; therefore, all specimens should be handled with the appropriate precautions.
8. After use, tubes and swabs must be disposed according to laboratory regulations for infectious waste.
9. The product must be used directly, it must not be subjected to any additional chemical or physical sterilization or micro-cidal or micro-static processes prior to use, as this will compromise the performance and function of the product.
10. If the intention is to use any part of a Copan Viral Transystem product with a third party test kit or assay then the user or manufacturer of such third party test kits or assays must verify acceptability of the Copan product or independently validate and qualify the use of the Venturi Transystem with said test kit or assay.
11.  Copan Viral Transystem is for single use only; reuse may cause a risk of infection and/or inaccurate results.
12. The performance is guaranteed only by the combination of the swab + transport tube hereby provided.
13. Do not pre-wet the collection swab into the transport medium prior to the execution of the collection.
14. Copan Viral Transystem medium does not contain antibiotics. Patient specimen that may contain a high load of bacterial contaminants may require additional antibiotics into the refeed medium.

PRODUCT DETERIORATION

The contents of unopened or undamaged units are guaranteed sterile. Do not use if they show evidence of damage, dehydration or contamination. Do not use if expiration date has passed.

INSTRUCTIONS FOR USE

Sample Transfer, Transport and Storage

- a) Peel open the Viral Transystem sterile pouch at the point marked "PEEL HERE."
- b) Remove cap from transport tube.
- c) Remove applicator swab holding from the plastic handle and collect specimen. During specimen collection, the applicator tip should only touch the area where the infection is suspected to minimize potential contamination.
- d) Place applicator swab in transport tube and close tightly.
- e) Record patient's name and information on tube label.
- f) Send specimen to the laboratory for immediate analysis.

Use in the laboratory

Copan Viral Transystem specimens should be processed for virology culture using recommended cells line and laboratory techniques that will depend on the specimen type and the organism under investigation. For recommended shell vials and techniques for the isolation and identification of viral strains from clinical swab specimens, refer to published virology manuals and guidelines. Culture investigations of swab specimens for the presence of viruses routinely involve the use of cells culture in shell vials. The procedure for inoculation of Copan Viral Transystem specimens onto shell vials is described below.

1. Note: Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens.
2. Squeeze using fingers and for 5 seconds the Copan Viral Transystem tube containing the swab specimen to release the specimen from the swab tip
3. Open the tube cap and remove the swab applicator.
4. Using a pair of scissor, cut the plastic tube where restricted. Squeeze the bottom of the tube with your fingers and, using a sterile pipette, move 200 µl volumes of the suspension into the shell vial or into a plate of a 24 wells plate, and proceed as per laboratory internal procedure.

NOTE: Patient specimen that may contain a high load of bacterial contaminants, may require additional antibiotics into the refeed medium.

QUALITY CONTROL

Copan Viral Transystem medium and applicators are tested to ensure they are non-toxic to cell lines used for viral cultivation. Copan Viral Transystem transport medium is tested for pH stability. Copan Viral Transystem is quality control tested before release for ability to maintain viable Coxsackie B1 ATCC VR28 virus. Procedures for quality control of virus transport medium and virus culture media are described in a number of publications by the American Society for Microbiology.

PERFORMANCE CHARACTERISTICS

Viability studies were performed using Copan Viral Transystem with a variety of viruses. Swabs accompanying each transport system were directly inoculated with 100 µL of virus suspension of known viral load. Swabs were then placed in their respective transport medium tubes and were held for 0, 6, and 24 hours at both 4°C and at room temperature (20-24°C). At the appropriate time interval, the tube was squeezed for at least 5 seconds using fingers, the swab was removed from its transport medium tube. Transport tube was cut and squeezed, then 200µL aliquot of the liquid medium was inoculated into shell vials. All cultures were processed by standard laboratory culture technique and tested by FA/DFSA. Organism viability was determined by fluorescing foci detection.

Viruses evaluated were: *Coxsackie Virus B5*, *Herpes Simplex Virus Type 1*, *Herpes Simplex Virus Type 2*, *Influenza A* and *Respiratory Syncytial Virus*. The results for the viruses tested using Copan Viral Transystem are shown in the table below.

Viral inoculated swabs samples	Viability at Room Temperature			Viability at Refrigerated Temperature		
	Holding Time	Incubation time of the cell culture before fixing, staining and reading (hours)	Detection of Foci of infected cells ¹	Holding Time	Incubation time of the cell culture before fixing, staining and reading (hours)	Detection of Foci of infected cells ¹
Coxsackie Virus B4	0 hr	24	TNTC	0 hr	24	TNTC
	6 hr	24	POS	6 hr	24	POS
	24 hr	24	POS	24 hr	24	POS
Herpes Simplex Virus Type 1	0 hr	24	229	0 hr	24	229
	6 hr	24	POS	6 hr	24	POS
	24 hr	24	POS	24 hr	24	POS
Herpes Simplex Virus Type 2	0 hr	24	256	0 hr	24	256
	6 hr	24	POS	6 hr	24	POS
	24 hr	24	POS	24 hr	24	POS
Influenza A	0 hr	24	264	0 hr	24	264
	6 hr	24	POS	6 hr	24	POS
	24 hr	24	POS	24 hr	24	POS
Varicella Zoster Virus	0 hr	24	TNTC ¹	0 hr	24	TNTC ¹
	6 hr	24	POS	6 hr	24	POS
	24 hr	24	POS	24 hr	24	POS
Respiratory Syncytial Virus	0 hr	24	285	0 hr	24	285
	6 hr	24	POS	6 hr	24	POS
	24 hr	24	POS	24 hr	24	POS

¹ TNTC = Too numerous to count, the complete monolayer is fluorescing.

POS = POSITIVE, triplicate tests performed on 200 µl aliquots of swab elution solution at each time point show positivity.

Copan Viral Transystem was able to maintain the viability of the following organisms for 24 hours at both room temperature (20-24°C) and in the refrigerator (2-8°C) under the test conditions described above: *Coxsackie virus B4*, *Herpes Simplex Virus Type 2*, *Influenza A*, *Respiratory Syncytial Virus*, *Varicella Zoster Virus*.

The refrigerated transport gave always better results than room temperature transport.

Performance testing with Copan Viral Transystem was conducted using pure laboratory virus strains spiked into the Copan Viral Transystem and not using clinical specimens. Proper specimen collection from the patient is extremely critical for the successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals. Specimens should be collected as soon as possible after the clinical onset of disease. Highest viral titers are present during the acute phase of the illness.

RESULTS

The amount of viable virus in a clinical specimen declines at a rate that varies depending on temperature and time in transit, the amount of virus originally present, the amount and type of protein in the specimen, the pH, bacterial flora¹. A specimen for virus isolation should be submitted to the laboratory as quickly as possible after collection. Copan Viral Transystem is able to maintain the viability of viruses for up to 24 hours.

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REFERENCES

- Lynne S. Garcia Clinical Microbiology Procedures Handbook, Third Edition 2010 ASM Press 3.2.1
- Steven Specter, Richard L. Hodinka, Stephen A. Young Clinical Virology Manual Third edition, 384-409
- Fleming D. Biological Safety: Principles and Practices. January 2000. ASM, Washington DC.
- Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA. <http://www.cdc.gov/od/ohs/symp5/jyrtext.htm>.
- Richardson JH. Biosafety in Microbiological and Biomedical Laboratories. December 1994. Diane Publishing Company.
- Hansen DJ. Healthcare, Laboratories and Biosafety. Vol 2, 1992. CRC Press.
- Gleaves C. A., R. L. Hodinka, S. L. G. Johnston and E. M. Swierkosz, Cumitech 15A. Laboratory Diagnosis of Viral Infections, p. 7. American Society for Microbiology, Washington D.C., 1994
- Miller, M. J., and A.L. Warford. Preparation of specimens for inoculation of cell cultures, p. 8.3.1 – 8.3.8. In H.D. Isenberg (ed.), Clinical Microbiology Procedures Handbook. American Society for Microbiology, Washington, D.C., 1992.
- Balows A., W.J. Hausler, Jr., K.L. Herrmann, H.D. Isenberg, H.J. Shadomy. Manual of Clinical Microbiology. Fifth Edition. American Society for Microbiology, Washington, DC, 1991.
- Wiedbrauk, D. L., and S. L. G. Johnston. Manual of Clinical Virology, p.11. Raven Press, New York, 1993.
- Johnson F. B., Transport of Viral Specimens, p. 120 – 131. Clinical Microbiology Reviews, Vol. 3, No. 2, April 1990.

Tabella dei Simboli / Table of Symbols / Tableau des pictogrammes / Tabelle der verwendeten Symbole / Tabla de símbolos / Tabela de Símbolos

Simbolo / Symbol / Pictogramme / Símbolo	Significato / Meaning / Signification / Bedeutung / Significado
	Fabbricante / Manufacturer / Fabricant / Hersteller / Fabricante
IVD	Dispositivo diagnostico in vitro / In vitro diagnostic device / Dispositif de diagnostic in vitro / In-vitro-Diagnostikum / Dispositivo de diagnóstico in vitro
CE 0123	Numero di identificazione dell'organismo notificato/ Identification number of notified body/ Identification de l'organisme notifié/ Identifizierung der benannten Stelle/ Identificación del organismo notificado/ Identificação do organismo notificado
STERILE R	Sterilizzato usando radiazioni ionizzanti / Sterilized using ionizing radiation / Stérilisé à l'aide de rayonnements ionisants / Sterilisation durch Bestrahlung / Esterilizado por radiaciones ionizantes / Esterilizado por radiação ionizante
	Non riutilizzare / Do not reuse / Ne pas réutiliser / Nicht zur Wiederverwendung / No reutilizar / Não reutilizar
REF	Numero di catalogo / Catalogue number / Numéro de catalogue / Bestellnummer / Número de catálogo
	Limiti di temperatura / Temperature limits / Limites de température / Temperaturbegrenzung / Límites de temperatura / Límites de temperatura
	Utilizzare entro / Use before / À utiliser avant / Verwendbar bis / Usar antes de / Usar antes de
	Consultare le istruzioni per l'uso / Consult the instructions for use / Consulter les instructions d'utilisation / Gebrauchsanweisung beachten / Consultar las instrucciones de uso / Consultar as instruções de utilização
	Strappare per aprire/ Peel/ Décoller/ Abziehen/ Desprender/ Destacável
LOT	Codice del lotto (partita) / Batch code (lot) / Code du lot (lot) / Chargenbezeichnung / Código del lote (Lote) / Código do lote (remessa)
	Contenuto sufficiente per <n> test / Contents sufficient for <n> tests / Contenu suffisant pour <n> tests / Inhalt ausreichend für <n> Prüfungen / Contenido suficiente para <n> pruebas / Conteúdo suficiente para <n> testes

Copan Italia SpA
Via Perotti, 10
25125 Brescia, Italy

Copan Italia Spa
Via Perotti 10
25125 Brescia Italy
Tel: +39 030 2687211
Fax: +39 030 2687250
E-mail: info@copangroup.com
Website: www.copangroup.com

Innovating Together™

North American Distributor:
Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562 USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com

This page is intentionally left blank