LBM[®] Selenite

Instructions for Use

Copan SELENITE Instruction for Use & User Guide

See the glossary of symbols at the end of the instruction for use

INTENDED USE

Copan **Selenite Broth** is a selective enrichment broth for clinical samples used for the culturing and the subsequent isolation of enteric pathogens including *Salmonella spp.* and *Shigella spp.*

SUMMARY AND PRINCIPLES

Diarrhoea can be defined as an increased frequency of bowel movements with at least three discharges a day accompanied by a decreased stool consistency. It is caused by a wide range of pathogens (bacteria, viruses and parasites) through different pathogenic mechanisms. *Salmonella spp.* and *Shigella spp.* infections come from the multiplication of intestinal bacteria. Most *Salmonella spp.* infections in humans are related to the consumption of animal based food, water contaminated by animals or interpersonal contact. Gastroenteritis is the most common infection that occurs. A small amount of *Salmonella spp.* may be present in the faeces of asymptomatic carriers. *Shigella spp.* infections have a wide variety of symptoms that range from mild diarrhoea to dysentery. Shigellas are transmitted directly from person to person, environmental multiplication is insignificant and asymptomatic infections are rare⁽²⁾. The aetiological diagnosis of *Salmonella spp.* and *Shigella spp.* diarrhoea is done by seeding the fecal sample directly onto a moderately selective solid culture medium such as Hektoen enteric agar/XLD (Xylose lysine deoxycholate), a more selective medium such as SS (Salmonella Shigella Agar)^(1,2) or by passing from an enrichment phase in a liquid medium such as Selenite broth^(1,2). Using an enrichment broth can increase the recovery of *Salmonella spp.* and *Shigella spp.* by 10% and is particularly recommended in populations such as service personnel and nutritionists⁽⁵⁾.

PRODUCT DESCRIPTION

Catalogue number.	Product description	Package dimensions
475CE.A	Polypropylene 12x80 mm skirted screw-cap tube with internal conical shape filled with 2 ml of Selenite broth.	50 units in each sales pack 6x50 units per box
475CE	Polypropylene 12x80 mm skirted screw-cap tube with internal conical shape filled with 2 ml of Selenite broth.	50 units per sales pack 6x50 units per case
4E021N.A	Polypropylene 12x80 mm round bottom tube with screw-cap filled with 2 ml of Selenite broth.	50 units per sales pack 6x50 units per case

FORMULA

Peptone
Lactose
Sodium phosphate salt
Sodium selenite salt

REQUIRED MATERIALS BUT NOT INCLUDED

Suitable materials for the isolation and culturing of bacteria; material for patient sample collection (different types of collection and transport systems, with Amies transport medium i.e. ESwab[®]/FecalSwab[™]).

Refer to the laboratory reference manuals for the recommended culture techniques and identification procedures i.e.

PRODUCT STORAGE

This product is ready for use and requires no further preparation. It must be stored in its original container at a temperature between 5° and 25°C until the time of use. Do not overheat. Do not incubate or freeze prior to use. Improper storage can cause a loss of effectiveness. Do not use after the expiration date, which is clearly printed on the outside of the box and on the test tube label.

LIMITATIONS

- The sample collecting conditions, times and volume for the culture are important variables for obtaining reliable culture results. Follow the recommended sample collection guidelines.
- The performance tests on the Copan Selenite broth were performed using ATCC[®] strains and not using human clinical samples.
- An adequate collection of samples from the patient is an extremely critical factor for the successful isolation and identification of infectious organisms. Consult the published reference manuals for specific guidance on the sample collecting procedures. Samples should be collected as soon as possible after clinical signs of the disease appear. Higher bacterial concentrations are present during the acute phase of the disease⁽⁷⁾.
- Selenite Broth is a selective and enrichment broth for the isolation of samples of *Salmonella spp.* and *Shigella spp.* The product was not designed for the enrichment of other genera and species of bacteria.
- Fecal samples require special attention. They must be taken in the early stages of the disease and must be cultured within 2 hours after collection. Under optimal conditions, the sample should be transferred into the transport medium immediately after collection (products such as the Copan Fecal Swab and Copan Cary Blair are suggested sample transport systems) in order to keep the organisms viable. Samples should never be refrigerated⁽⁷⁾.
- There is no dulcitol (melampirite or polyalcohol) in the Copan Selenite Broth therefore the growth of coliforms could be suppressed during the first 8-12 hours of incubation; longer incubation times may result in the excessive growth of mixed flora⁽⁸⁾.

WARNINGS AND PRECAUTIONS FOR USE

- Observe approved biohazard precautions and aseptic techniques.
- The product is only to be used by adequately trained and qualified personnel.
- All samples and the materials used to process them must be considered potentially infective and must be handled in such a way as to prevent infection of the laboratory staff.
- Observe the other applicable biosafety recommendations.

5. For *in vitro* diagnostic use. This product must not be used for therapeutic or prophylactic purposes.
6. Carefully read and follow the instructions.
7. This product is intended for single use; if reused, there is the risk of obtaining misleading results.
8. Not suitable for any other application different from its intended use.
9. Work under a biosafe laminar flow hood. Wear latex gloves and other protection commensurate with universal precautions when handling clinical specimens.
10. The use of this product in association with any diagnostic tests or any diagnostic instrumentation should be previously evaluated by the user.
11. Do not use the product if it is visibly damaged.
12. The liquid contained in the individual tubes is TOXIC. Do not swallow or inhale and avoid contact with the skin and eyes. Take appropriate protective measures.
13. All samples are to be considered potentially infectious and must be handled with appropriate precautions.
14. Consult the operating instructions supplied with the device or available in electronic format, and which can be identified by the e-IFU indicator on the packaging label.
15. The product is only to be used by properly trained and qualified personnel.
- 16.

Warning: Contains Sodium Selenite - Harmful if swallowed

Wash hands thoroughly after handling - Do not eat, drink or smoke when using this product - IF SWALLOWED: if you feel unwell call a POISON CENTER/ doctor - Rinse mouth - Dispose of contents/container in accordance with local regulation.

PRODUCT DETERIORATION

Do not use Copan Selenite Broth if: 1. the product has visible signs of damage or contamination; 2. there is evidence of liquid leaking from the test tube; 3. it has passed the expiration date; 4. there are other visible signs of deterioration.

INSTRUCTIONS FOR USE

Use in the laboratory

Selenite Broth is used for the selective enrichment of *Salmonella spp* and *Shigella spp* present in fecal or rectal samples. The Selenite Broth comes with a special cap that has an adapter for attaching the shaft of the swab used during sampling (sampling swabs with a marked break point on the shaft are recommended for use with the Selenite Broth tubes. For further information please contact Copan Italia S.p.A.).

If a rectal or fecal swab is obtained using the FecalSwab™ collection and transport system, this swab can be transferred into the Selenite Broth as indicated in the procedure below.

1. Unscrew the cap of the Selenite Broth tube.
2. Add the sample as described below:

Rectal or fecal swab: (homogenize the sampling system with the sample by vortexing the test tube containing the sample for 10 seconds at 2000-2500 rpm).

- For a rectal or fecal swab obtained using FecalSwab™: open the tube and transfer the swab into the Selenite Broth using sterile tweezers; alternatively, transfer approximately 30µl of culture medium in the FecalSwab™ system using a micropipette; or transfer the shaft attached to the cap of the FecalSwab™ system directly into the Selenite Broth, then close the FecalSwab™ tube with the cap from the Selenite tube.
- For a generic rectal or fecal swab or a dry swab (without transport medium): break the swab at the marked break point inside the Selenite Broth and close the cap, or vigorously shake the swab without a break point in the tube for at least 10 seconds.
- Close the inoculated Selenite Broth tube.
- Recap the tube/container with the original sample.

Fecal sample:

- Add a portion of the fecal sample using a swab or a special scoop.
NOTE: make sure that the swab used to transfer the feces is completely covered with fecal material.
 - Transfer a portion of the feces into the Selenite Broth. Follow the guidelines used in your laboratory for the volumetric ratio.
 - Vigorously shake the swab or scoop in the Selenite Broth for at least 10 seconds, in order to release the collected material.
 - Close the Selenite Broth tube.
3. Homogenize the tube by vortexing it for 10 seconds at 2000-2500 rpm;
 4. Incubate the Selenite Broth tubes at 35°C ± 2°C for 18-24 hours;
 5. After 18-24 hours look for any turbidity in the broth. If negative, incubate for another 18-24 hours;
 6. Inoculate a plate of selective and/or differential growth media with 1 to 10 µl of enriched broth.

The enriched sample in Selenite Broth can be seeded with automated sample management systems such as WASP®/WASPLab®. For seeding, follow the procedure for seeding broths given in the WASP® automation manual.

QUALITY CONTROL

ENRICHMENT VERIFICATION

- Prepare a 0.5 Mac Farland inoculum in PBS from a fresh culture of *Salmonella enterica* ATCC® 14028.
- Prepare a 10⁻⁵ in series dilution from the 0.5 Mac Farland in order to obtain a plate count of between 30 and 300 CFU/100 µL.
- Open the tube of Selenite Broth under a laminar flow hood.
- Using a micropipette, inoculate 200 µL of the suspension directly in the 2 ml tube of Selenite Broth (or 500 µL in a 5 ml tube).
- Recap the tube.
- Vortex for 10 seconds at 2000-2500 rpm.
- Immediately re-open the tube and streak an XLD agar plate with 100 µL or following the internal laboratory procedure.
NOTE: The plate must be marked as time zero for the test.
- Incubate the tube at 35°C ± 2°C for 18-24 hours.

- After incubation, streak an XLD agar plate with 100 µL of Selenite Broth or following the internal laboratory procedure using another type of culture medium suitable for isolating *Salmonella* spp.
- Incubate the agar plate as required in the standard laboratory procedure, at 35°C ± 2°C for 18-24 hours.

The plate count at time zero must be 30-300 CFU/plate and confluent growth must be visible after 18-24 hours.

INHIBITION TEST:

- Prepare a 0.5 Mac Farland inoculum in PBS from a fresh culture of *E. coli* ATCC® 25922.
- Open the tube of Selenite Broth under a laminar flow hood.
- Using a micropipette, inoculate 200 µL of the inoculum directly in the 2 ml tube of Selenite Broth (or 500 µL in a 5 mL tube).
- Recap the tube.
- Vortex for 10 seconds at 2000-2500 rpm.
- Immediately re-open the tube and streak a plate of McConkey agar with 100 µL or following the internal laboratory procedure using another type of culture medium suitable for isolating *E. coli*. The plate must be marked as time zero for the test.
- Incubate at 35°C ± 2°C for 18-24 hours.
- After incubation, streak a McConkey agar plate with 100 µL of Selenite Broth or following the internal laboratory procedure using another type of culture medium suitable for isolating *E. coli*.
- Incubate the agar plate as required in the standard laboratory procedure, at 35°C ± 2°C for 18-24 hours.

The plate must show confluent growth at time zero and partial or complete inhibition must be visible after 18-24 hours.

Strains	Results CFU/plate at time zero Average of 5 plates	Results CFU/plate after 24 hours Average of 5 plates
<i>Salmonella enterica</i> Serovar Typhimurium ATCC® 14028	172	CONFLUENT GROWTH
<i>E. coli</i> ATCC® 25922	CONFLUENT GROWTH	PARTIAL TO TOTAL INHIBITION
<i>Shigella sonnei</i> ATCC® 9290	202	CONFLUENT GROWTH

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BIBLIOGRAPHY

1. J.H. Jorgensen, M.A. Pfaller, K.C.Carroll, G.Funke, M.L.Landry, S.S.Richter,., 2015. Manual of Clinical Microbiology.11 th ed. ASM, Washington, D.C.
2. UK Standards for Microbiology Investigations. of fecal specimens for Enteric Pathogens. Bacteriology | B 30 | Issue no: 8.1 | Issue date: 24.04.14.
3. Clinical and Laboratory Standards Institute (CLSI). 2004. Quality Control for commercially Prepared Microbiological Culture Media. Approved Standard – Third Edition M22-A3
4. Miller, J. M. 1999. A Guide to Specimen Management in Clinical Microbiology, 2nd ed. ASM, Washington, DC.
5. Garcia L.S. H. D., 2010. Clinical Microbiology Procedures Handbook, 3rd ed. ASM, Washington, DC.
6. Atlas, R.N.1993 Handbook of Microbiological Media p. 147-153. CRC Press, Boca Raton, FL
7. P.H. Gilligan,J.M.Janda, M.A. Karmali, J.M. Miller. Laboratory Diagnosis of Bacterial Diarrhea. Cumitech 12A, American Society for Microbiology, Washington DC 1992.
8. E. W. Koneman,S.D.allen, V.R. Dowell Jr., W.M. Janda, H.M. Sommers, W. C. Winn Jr. The Enterobacteriaceae in: Diagnostic Microbiology. Third Ed. J.B. Lippincott Company 1988.

INDEX OF SYMBOLS

Symbol / Simbolo / Símbolo / Symbole / Symbole / Símbolos / Symbol / Symbol / Símbol / Sembol	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Betyder / Betydning / Betydelse / Anlami
	Manufacturer / Fabbicante / Fabricante / Hersteller / Fabricant / Fabricante / Producent / Produsent / Tillverkare / Üretici
	In vitro diagnostic device / Dispositivo Diagnostico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro/ Dispositivo de diagnóstico in vitro / diagnostisk medicinsk anordning in vitro / In vitro diagnostisk medisinsk utstyr / Medicinsk anordning / För in vitro-diagnostik
	CE marking / Marchio CE / Marca CE / CE-Kennzeichnung / Marquage CE/ Marcação CE / CE-mærke / CE-merke / CE-märkning / CE işareti
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Brug det ikke igen / Må ikke gjenbrukes / Ateranvänd inte / Tekrar kullanmayın
	Catalogue number / Numero di catalogo / Número de catálogo / Bestellnummer / Référence du catalogue / Número do catálogo / Katalognummer / Katalognr / Katalognummer / Katalog numarası
	Temperature limitation / Limiti di temperatura / Límites de temperatura / Temperature Begrenzung / Límites de temperatura / Límites de temperatura / Temperaturgrænser / Temperaturgrenser / Temperaturgrænser / Sicaklık limitleri
	Use by / Utilizzare entro / Fecha de caducidad / Verwendbar bis / Utiliser jusque / Prazo de validade / Anvendes for / Må brukes innen / Sista förbrukningsdag / Son kullanma tarihi
 eIFU Indicator	Consult the operating instructions supplied with the device or available in electronic format, and which can be identified by the e-IFU indicator on the packaging label / Consultare le istruzioni per l'uso fornite con il dispositivo oppure disponibili in formato elettronico ed identificate dall'e-IFU indicator sull'etichetta imballo / Consultar las instrucciones de uso suministradas con el dispositivo o disponibles en formato electrónico e identificadas por el indicador e-IFU de la etiqueta del embalaje / Konsultieren Sie die Gebrauchsanweisung, die mit dem Gerät geliefert wird oder in elektronischem Format vorliegt und durch den e-IFU-Indikator auf dem Verpackungsetikett gekennzeichnet ist / Voir le mode d'emploi fourni avec l'appareil ou disponible au format électronique et identifiable grâce à l'indicateur e-IFU sur l'étiquette de l'emballage / Ver as instruções de utilização fornecidas com o dispositivo ou disponíveis em formato eletrônico e identificadas pelo indicador e-IFU na etiqueta da embalagem / Se brugsvejledningen, der følger med enheden eller er tilgængelig i elektronisk format og identificeret ved e-IFU-indikatoren på emballagemærkatet/ Se bruksanvisningen eller e-bruksanvisningen hvis «eIFU-indikator» er til stede / Se bruksanvisningen som bifogas apparaten eller finns tillgänglig i elektronisk form via IFU-koden på förpackningen / Cihaz ile tedarik edilen veya elektronik şekilde ve ambalaj etiketinde e-IFU kodu ile tespit edilen kullanım talimatlarına bakınız
	Batch code (Lot) / Codice del lotto (partita) / Código de lote (Lote) / Chargencode (Chagenbezeichnung) / Code de lot (Lot) / Código do lote (Lote) / Serienummer (parti) / Lot-nummer (parti) / Partnummer (lotto) / Lot (seri) kodu
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> test / Contenido suficiente para <n> pruebas / Ausreichend für <n> Tests / Contenu suffisant pour <n> tests / Contém o suficiente para <n> testes / Indhold tilstrækkeligt til <n> prøver / Innhold tilstrekkelig for <n> prøver / Innehåller tillräckligt för <n> tester / <n> test için yeterli içerik
	Not-sterile / Non sterile / Non stérile / Nicht steril / No estéril / Não estéril / Ikke steril / Ikke steril / Ikke-steril / Steril değildir
	Caution / Attenzione / Precaución / Vorsicht / Attention / Cuidado / Forsigtig / Forsiktig / Laktag försiktighet / Dikkat

Copan Italia S.p.A.
Via F. Perotti, 10
25125 Brescia, Italy
Tel +39 030 2687211
Fax +39 030 2687250

Email: info@copangroup.com
Website: www.copangroup.com

North American Distributor:

Copan Diagnostics Inc.
26055 Jefferson Avenue
Murrieta, CA 92562, USA
Tel: 951-696-6957
Fax: 951-600-1832

E-mail: customerservice@copanusa.net
Website: www.copanusa.com