UriSponge™

Instructions for Use

IVD CE 0123

English	3
Italiano	6
Español	9
Deutsch	13
Français	16
Português	19
БЪЛГАРСКИ	23
Česky	26
Dansk	30
ΕΛΛΗΝΙΚΑ	33
Eesti Keel	37
Hrvatski	40
Latviešu	43
Lietuvių K.	46
Magyar	50
Nederlands	53
Norsk	57
Polski	60
Română	63
Slovenčina	66
Slovenščina	69
Suomi	73
Svenska	76
Türkçe	79

Copan UriSponge™ - Urine Collection, Transport and Preservation System
Instructions for use.

INTENDED USE

Urine transport and preservation system

Copan UriSponge™ - Urine Collection, Transport and Preservation System is intended for the transport and preservation of urine specimens, transferred from their initial container, from the collection site to the testing laboratory. In the laboratory, UriSponge™ specimens are processed using standard clinical laboratory operating procedures for the cultivation of uropathogenic bacteria and yeasts.

SUMMARY AND PRINCIPLES

Quantification of the bacterial load of uncontaminated intermediate urination is a widely used system to evaluate urinary tract infections¹. Maintaining the bacterial load in urine specimens depends on numerous factors, including the type and concentration of microorganisms, transport duration and storage temperature². The UriSponge™ device contains preservative substances (urine preservation formula) on the applicator sponge in a conical bottom test tube with screw cap that preserve the viability of clinically important organisms during transport for up to 48 hours at 2-25°C.

REAGENTS

Copan UriSponge™ contains the following preservatives (urine preservation formula) on the applicator sponge:

- Sodium propionate
- Potassium sorbate

PRODUCT DESCRIPTION

The Copan UriSponge™ collection, transport and preservation system is supplied ready for use. The product descriptions and packaging configurations are listed in Table 1.

The Copan UriSponge™ collection, transport and preservation system is sterilized by ionising radiation.

The codes with suffix "A" can be processed with automated systems. Refer to the manufacturer's instructions on how to use UriSponge™ with the automated system.

Figure 1: UriSponge™ device


Figure 2: UriSponge™ components

Table 1. UriSponge™ description and packaging

Catalogue number	Product description	Packaging
802CE.A	Urine transport and preservation system, 12x80 test tube	50 devices per sales package 50 x 6 devices per box

PRODUCT STORAGE

Store in its original package at a temperature between 2 and 25°C until the time of use. Do not freeze before use.

REQUIRED MATERIALS BUT NOT INCLUDED

Sterile urine collection container and appropriate materials for isolation and culture of uropathogenic bacteria and yeasts, e.g. culture media and incubation systems are not provided. For the recommended procedures for culture and identification of uropathogenic bacteria and yeasts of clinical specimens, refer to the reference laboratory manuals^{1,4}.

SPECIMEN STORAGE AND TRANSPORT

UriSponge™ specimens should immediately be transferred to the laboratory, preferably within 2 hours of collection². If immediate delivery or analysis is delayed, the specimens collected with UriSponge™ must be refrigerated at 2-8°C or stored at room temperature (20-25°C) and analysed within 48 hours.

LIMITATIONS

1. The conditions, times and volume of the specimen collected for culture are significant variables for obtaining reliable culture results. Follow the recommended guidelines for specimen collection^{1,2,3,4}.
2. The microbial load in the urine of a given patient may be affected by the time of collection and the amount of liquids taken. Symptomatic patients may have counts of less than 10E+05 microorganisms/ml if specimens are collected late in the afternoon or in the presence of diuresis¹⁰.
3. The preservatives contained in the sponges do not inactivate antibiotics.
4. Not suitable for recovering microorganisms that are difficult to handle, such as anaerobes, viruses, chlamydias, mycoplasmas, ureaplasma, trichomonas, in culture.
5. Not suitable for collection of specimens for chemical-physical, microscopic and macroscopic examinations.
6. Not validated for use with molecular methods.

WARNINGS and PRECAUTIONS

1. For in vitro diagnostic use only.
2. For professional use only. Use UriSponge™ as indicated in the instructions for use.
3. Do not use if the package or test tube is open, if the product shows signs of damage, deterioration or contamination or if it has passed the expiry date. Do not touch the absorbent sponges.
4. Do not reuse. Do not resterilize. Do not repackage.
5. The product is sterile; do not use the product if the sterile barrier is damaged.
6. The use of this device in association with diagnostic instruments and/or kits must be validated prior to use.
7. Wear gloves and other protection devices following the general precautions for handling of clinical specimens. Observe the CDC biosafety recommendations. After use, dispose of the test tubes according to the laboratory regulations for hazardous waste^{7,8,9}. Before transport, make sure that the UriSponge™ screw cap is tightly closed.
8. Check the version of the instructions for use. The correct version is the one provided with the device or available in electronic format and identified by the e-IFU indicator on the package label.

INSTRUCTIONS FOR USE - COLLECTION

Specimen collection and transport with UriSponge™

Proper collection of patient specimens is an essential element to ensure the isolation and identification of infectious organisms. Specimens obtained for evaluation of urinary tract pathogens must be collected and handled according to published manuals and guidelines^{1,2,3,4}. Collect a urine specimen from intermediate urination into a sterile container.

Open the UriSponge™ test tube and, holding it by the cap, immerse the sponge applicator in the urine specimen. Keep the sponges immersed for not less than 5 seconds. The polyurethane sponges are extremely hydrophilic and spontaneously absorb the urine.


Remove the sponge applicator from the urine specimen and re-insert it into the UriSponge™ test tube. See Figure 3. Do not add urine to the test tube other than that spontaneously absorbed by the sponge.

Screw on the cap to securely close the container.

Label and transport the specimens collected with UriSponge™ in accordance with institutional, local, state and federal requirements⁶.

Immediately transfer the specimens collected with UriSponge™ to the laboratory, preferably within 2 hours of collection². If immediate delivery or analysis is delayed, the specimens collected with UriSponge™ must be refrigerated at 2-8°C or stored at room temperature (20-25°C) and analysed within 48 hours.

Figure 3: Transfer of urine specimen to the UriSponge™ device.

INSTRUCTIONS FOR USE - SPECIMEN PROCESSING IN THE LABORATORY

Analysis of specimens collected with UriSponge™ for laboratory culture

Specimens collected with UriSponge™ must be analysed immediately upon receipt and within 48 hours of collection.

The UriSponge™ test tubes can be centrifuged to extract the urine from the sponge. Centrifuge at 104g in case of a 12x80 test tube for 3 seconds in a rotor centrifuge. The UriSponge™ test tubes can be re-centrifuged as needed. Alternatively, the specimen can be manually extracted from the sponge by shaking the device three times with rapid, decisive downward wrist movements.

Using an aseptic technique, unscrew the cap and remove the sponge applicator. Gently mix the contents of the test tube.

Follow the laboratory internal standard operating procedures (SOPs) to inoculate the urine specimen of the test tube on a culture medium or refer to published microbiology manuals and guidelines^{1,4}.

QUALITY CONTROL

UriSponge™ has been verified to preserve viability after storage at 2-8°C (refrigeration) and at 20-25°C (room temperature) for up to 48 hours using the organisms recommended by the Clinical and Laboratory Standards Institute (CLSI) M40-A2⁵.

The procedures for quality control of microbiological specimen transport devices must be conducted using the analysis methods described in CLSI M40-A2⁵. If abnormal quality control results are observed, patient results must not be reported.

PERFORMANCE CHARACTERISTICS

Tests were conducted using the organisms recommended by the Clinical and Laboratory Standards Institute (CLSI) M40-A2⁵ and other ATCC® strains in artificial urine

Urisponge™ was saturated according to the instructions for use.

Urisponge™ was verified to preserve viability after storage at 2-8°C (refrigeration) and at 20-25°C (room temperature) for up to 48 hours.

Table 2. Urisponge™ performance in preserving the viability of urinary tract microorganisms

CLSI M40 strains:

Organism	Holding Temperature	Time zero (CFU/mL)	Time 24 hours (CFU/mL)	Logarithmic difference of Time 24h	Time 48 hours (CFU/mL)	Logarithmic difference of Time 48h
<i>C. albicans</i> (ATCC® 24433)	2-6°C	8,00E+02	6,60E+02	-0,08	5,40E+02	-0,17
	20-24°C	8,00E+02	9,60E+02	0,08	1,77E+03	0,35
<i>E. coli</i> (ATCC® 25922)	2-6°C	1,17E+03	7,80E+02	-0,18	8,36E+02	-0,15
	20-24°C	1,17E+03	5,10E+02	-0,36	8,37E+02	-0,15
<i>E. faecalis</i> (ATCC® 29212)	2-6°C	1,44E+03	7,00E+02	-0,31	7,10E+02	-0,31
	20-24°C	1,44E+03	1,71E+03	0,07	1,90E+03	0,12
<i>P. aeruginosa</i> (ATCC® 27853)	2-6°C	6,90E+02	6,60E+02	-0,02	4,90E+02	-0,15
	20-24°C	6,90E+02	4,00E+02	-0,24	1,20E+02	-0,76
<i>P. mirabilis</i> (ATCC® 7002)	2-6°C	1,22E+03	1,04E+03	-0,07	5,90E+02	-0,32
	20-24°C	1,22E+03	9,40E+02	-0,11	6,93E+02	-0,25
<i>S. saprophyticus</i> (ATCC® 15305)	2-6°C	9,70E+02	9,90E+02	0,01	6,90E+02	-0,15
	20-24°C	9,70E+02	9,30E+02	-0,02	6,23E+02	-0,19

Additional ATCC® strains

Organism	Holding Temperature	Time zero (CFU/mL)	Time 24 hours (CFU/mL)	Logarithmic difference of Time 24h	Time 48 hours (CFU/mL)	Logarithmic difference of Time 48h
<i>C. freundii</i> (ATCC® 8090)	2-6°C	1,22E+03	2,60E+02	-0,67	1,60E+02	-0,88
	20-24°C	1,22E+03	8,55E+02	-0,15	6,70E+02	-0,26
<i>C. glabrata</i> (ATCC® 15126)	2-6°C	2,30E+03	1,41E+03	-0,21	1,10E+03	-0,32
	20-24°C	2,30E+03	1,87E+03	-0,09	2,40E+03	0,02
<i>E. cloacae</i> (ATCC® 13047)	2-6°C	1,54E+03	1,24E+03	-0,09	1,28E+03	-0,08
	20-24°C	1,54E+03	1,06E+03	-0,16	1,01E+03	-0,18
<i>M. morgani</i> (ATCC® 25829)	2-6°C	3,00E+03	1,88E+03	-0,20	2,15E+03	-0,14
	20-24°C	3,00E+03	2,15E+03	-0,15	3,13E+03	0,02
<i>K. pneumoniae</i> (ATCC® 700-603)	2-6°C	1,58E+03	7,80E+02	-0,31	6,20E+02	-0,41
	20-24°C	1,58E+03	9,50E+02	-0,22	8,20E+02	-0,28

Microorganism release study (centrifugation and manual release methods)

The ability of Urisponge™ to release suspended organisms in filtered synthetic urine at known concentrations by a centrifugation and manual method was compared with the concentration of the organisms in the initial inoculated synthetic urine. See Table 3.

Table 3. Performance of Urisponge™ in releasing a selected series of uropathogenic microorganisms using the centrifugation and manual release method

Organism	Mean CFU/ml Recovered: Inoculated urine	Mean CFU/ml Recovered: Time 0 hours Release by centrifugation	Mean CFU/ml Recovered: Time 0 hours Manual release	Release by centrifugation vs. Logarithmic reduction (-) or increase (+) in inoculated urine	Manual release vs. Logarithmic reduction (-) or increase (+) in inoculated urine
<i>C. albicans</i> (ATCC® 24433)	7.70E+02	7.70E+02	9.90E+02	0.00	0.11
<i>E. coli</i> (ATCC® 25922)	7.80E+02	1.26E+03	1.12E+03	0.21	0.16
<i>S. saprophyticus</i> (ATCC® 15305)	1.13E+03	9.90E+02	1.06E+03	-0.06	-0.03

NOTES FOR THE PROFESSIONAL USER

In the event that a serious incident occurs in relation to this device, it must be reported to the manufacturer (see the contacts at the end of the Instructions for Use) and the competent authority in the country where the user and/or patient is located.

REVISION HISTORY

Last Revision No.*	Date of issue	Changes made
02	07-2023	Update for new formulation of the salts contained in the sponge

* Should you need earlier revisions, contact Copan Customer Service.

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BIBLIOGRAPHY

1. Yvette S. Mccarter, Eileen M. Burd, Gerri S. Hall, And Marcus Zervos Laboratory Diagnosis of Urinary Tract Infections. Cumitech 2C, 2009. ASM, Washington DC.
2. Miller JM. A Guide to Specimen Management in Clinical Microbiology. Third Edition. American Society for Microbiology. Washington, DC. 2017.
3. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 2001. Urinalysis and Collection, Transportation, and Preservation of Urine Specimens, Approved Guideline. Second Edition. GP16-A2. Vol.21 No.19.
4. Lynne S. Garcia, Clinical Microbiology Procedures Handbook, Third Edition ASM, Washington DC 2010.
5. Clinical Laboratory Standards Institute CLSI (formerly National Committee for Clinical Laboratory Standards NCCLS). 2014. Quality Control of Microbiological Transport Systems; Approved Standard. Second Edition M40-A2 Vol. 34 No. 9.
6. 542CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
7. Richard J. The 1, 2, 3's of Biosafety Levels. Centers for Disease Control and Prevention, Atlanta, GA.
8. Biosafety in Microbiological and Biomedical Laboratories (BMBL) 6th Edition, U.S. Department of Health and Human Services Public Health Service Centers for Disease Control and Prevention National Institutes of Health, Revised June 2020
9. LaRocco MT et al. Effectiveness of Preanalytic Practices on Contamination and Diagnostic Accuracy of Urine Cultures: a Laboratory Medicine Best Practices Systematic Review and Meta-analysis. Clin Microbiol Rev. 2016 Jan;29(1):105-47. Doi: 10.1128/CMR.00030-15. PMID: 26598386; PMCID: PMC4771218.

TABLE OF SYMBOLS / TABELLA DEI SIMBOLI / TABLA DE SÍMBOLOS / SYMBOLTABELLE / TABLEAU DES SYMBOLES / TABELA DE SÍMBOLOS / ТАБЛИЦА СЪС СИМВОЛИ / TABULKA SYMBOLŮ / TABEL OVER SYMBOLER / ΠΙΝΑΚΑΣ ΣΥΜΒΟΛΩΝ / SÚMBOLITE TABEL / TABELICA SIMBOLA / SIMBOLU TABULA / SIMBOLIUR LENTELE / SZIMBÓLUMTÁBLÁZAT / TABEL MET SYMBOLEN / SYMBOLTABELL / TABELA SIMBOLA / TABELUL SIMBOLURILOR / TABULKA SYMBOLOV / TABELA SIMBOLOV / SYMBOLITAUULKKO / SYMBOLFÖRTECKNING / SEMBOLLER TABLOSU

Symbol / Simbolo / Símbolo / Symbol / Symbole / Símbolos / Символ / Symbol / Symbol / Σύμβολο / Tingmärk / Simbol / Simbols / Simbolis / Szimbólum / Symbol / Symbol / Symbol / Símbol / Symbol / Simbol / Simboli / Symbol / Sembol	Meaning / Significato / Signification / Bedeutung / Sens / Significado / Значение / Význam / Betydning / Σημασία / Tāhendus / Značenje / Nozīme / Reikšmė / Jelentés / Betekenis / Betydning / Znaczenie / Semnificație / Význam / Pomen / Merkitys / Betydelse / Anlam
	Manufacturer / Fabricante / Fabricante / Hersteller / Fabricant / Fabricante / Производител / Výrobce / Producent / Κατασκευαστής / Toolja / Proizvođač / Ražotājs / Gamintojas / Gyártó / Fabrikant / Produzent / Producător / Výrobca / Proizvajalec / Valmistaja / Tilverkare / Üretici
	In vitro diagnostic device / Dispositivo Diagnóstico in Vitro / Dispositivo de diagnóstico in vitro / In-vitro-Diagnostikum / Dispositif de diagnostic in vitro / Dispositivo de diagnóstico in vitro / Инвитро диагностическое изделие / Diagnostický prostředek pro in vitro diagnostiku / Udstyr til in vitro-diagnostik / Διαγνωστικό προϊόν in vitro / In vitro diagnostika seade / In vitro dijagnostički proizvod / In vitro diagnostikas ierīce / „In vitro“ diagnostikos prietaisais / In vitro diagnostikai eszköz / Hulpmiddel voor in-vitrodiagnostiek / In vitro-diagnostisk enhet / Urządzenie do diagnostyki in vitro / Dispozitiv de diagnostic in vitro / Diagnostická pomůcka in vitro / In vitro diagnostični pripomoček / Dijagnostički uređaj u epruveti / Enhet för in vitro-diagnostik / In vitro tani cihaz
	Identification number of notified body / Numero di identificazione dell'organismo notificato / Identificación del organismo notificado / Identifizierung der benannten Stelle / Identification de l'organisme notifié / Identificação do organismo notificado / Идентификационный номер на нотифицирующего орган / Identifikační číslo notifikovaného subjektu / Identifikationsnummer af bemyndiget organ / Αριθμός ταυτοποίησης κοινοποιημένου οργανισμού / Teavitatud asutuse identifitseerimisnumber / Identifikacijski broj obavještenog tijela / Plinvarotās iestādes identifikācijas numurs / Notifikuotosios įstaigos identifikavimo numeris / A bejelentéssel szerezett azonosító száma / Identificatienummer van de aangemelde instantie / Identifikasjonsnummeret til det tekniske kontrollorganet / Numer identyfikacyjny jednostki notyfikowanej / Numărul de identificare al organismului notificat / Identifikačné číslo notifikovanej osoby / Identifikacijska številka priglašenej organa / Ilmoitetun laitoksen tunnusnumero / Identifikationsnummer för anmält organ / Onaylı kuruluş tanımlama numarası
	Sterilized using irradiation / Sterilizzato usando radiazioni ionizzanti / Esterilizado usando radiaciones ionizantes / Sterilisiert mit ionisierenden Strahlungen / Stérilisé à l'aide de radiations ionisantes / Esterilizado por radiação ionizante / Стерилизовано чрез облъчване / Sterilizováno ozařováním / Steriliseret med stråling / Αποστειρωμένο με ακτινοβολία / Steriliseeritud, kasutades kiirust / Sterilizirano zračenjem / Sterilizēts, izmantojot apstarošanu / Sterilizuota švitināt / Sugárzással sterilizálva / Gesteriliseerd met bestraling / Steriliseret ved stråling / Sterylizowany promieniowaniem / Sterilizat prin iradiere / Sterilizované ožiarím / Sterilizirano z ionizirajočim sevanjem / Steriloitu säteilyttämällä / Steriliserad med stråling / Įšnįlama kulanlarak sterilize edilmıştır
	Sterile barrier system / Sistema di barriera sterile / Sistema de barrera estéril / Sterilbarrieresystem / Système de barrière stérile / Sistema de barreira estéril / Стерилна бариерна система / System sterilin bariéry / Sterilt barrieresystem / Σύστημα στείρου φραγμού / Ühekordne steriilne barjäärüsteem / Sustavi sterilne barijere / Sterilas barjeras sistema / Sterilaus barjero sistema / Steril gätrendszer / Steril barrieresysteem / Sterilt barrieresystem / System bariery sterylnej / Sistem de barieră steril / Systeeme sterilnej bariéry / Sterilni pregradni sistem / Steriili suojajärjestelmä / Sterilt barriersystem / Steril ambalaj sistemi
	Do not resterilize / Non risterilizzare / No reesterilizar / Nicht resterilisieren / Ne pas restériliser / Não reesterilizar / Не стерилизуйте повторно / Výrobek opakovaně nesterilizujte / Må ikke gensteriliseres / Μην επαναποστειρώνετε / Mitte resteriliseerida / Ne sterilizirajte ponovno / Nesterilizēt atkārtoti / Nesterilizuokite pakartotinai / Ne sterilizálja újra / Niet opnieuw steriliseren / Må ikke steriliseres på nytt / Nie sterylizować ponownie / A nu se resteriliza / Nesterilizujte znovu / Ne sterilizirajte ponovno / Ei saa steriloida uudelleen / Får ej omsteriliseras / Tekrar sterilize etmeyin
	Do not reuse / Non riutilizzare / No reutilizar / Nicht zur Wiederverwendung / Ne pas réutiliser / Não voltar a utilizar / Не используйте повторно / Výrobek nepoužívajte opakovaně / Må ikke genbruges / Μην επαναχρησιμοποιείτε / Ühekordne kasutus / Ne upotrebľajte ponovno / Nelietot atkārtoti / Nenaudoti pakartotinai / Ne használja újra / Niet hergebruiken / Skal ikke brukes flere ganger / Nie używać ponownie / A nu se reutiliza / Nepoužívajte opakovane / Ne uporabljajte ponovno / Ei saa käyttää uudelleen / Får ej återanvändas / Yeniden kullanılmayın

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Email: info@copangroup.com
 Website: www.copangroup.com

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Copan Italia S.p.A.