Product Focus

LBM® Broths

LBM® Broths

Liquid-Based Microbiology™ for every need





Among transport and processing media, the Copan LBM® family includes a full range of enrichment and selective broths for the most common aerobic and anaerobic bacteria. All LBM® broth tubes have been conceived with automation in mind and are designed explicitly to be loaded and processed by WASP®.



Durable

Convenient shelf-life and stability at room temperature (20-25°C)



Ready-to-use

No pretreatment and preparation needed



Compliant

All the LBM® broths family is (EU) 2017/746 (IVDR), ISO 13485-2016, ISO 9001-2015 compliant



Standardized

LBM® broths screw cap tube format fits most automated devices and laboratory equipment



Fields of application

Preanalytics made different



Gastrointestinal diseases



Antibiotic resistances



Women health (UTI and GBS)



Respiratory diseases



Hospital-acquired infections



FLOQSwabs®

LBM broths best ally

FLOQSwabs® offer variable sizes, diameters, breaking points and tip shapes to be used in plenty of applications. This made FLOQSwabs® a well-tolerated alternative to invasive, painful, and costly collection procedures^{7,8}.

Do you have a specific application in mind? Choose the right FLOQSwabs®!



WASP®

Walk-Away Specimen Processor™

Copan WASP® is a truly revolutionary instrument of specimen processing for Microbiology.

WASP® provides a comprehensive system encompassing all aspects of automated specimen processing, planting and streaking, Gram slide preparation, and enrichment broth inoculation.

Ordering information

Select the LBM® broth you need, in bulk packs or in combination with a swab (available only for specific products).

Cat N. Description Pack size Sample*

BHI (Brain Heart Infusion) Enrichment Broth™: BHI broth™ is a medium specifically intended for fastidious bacteria, including Streptococcus pneumoniae and Staphylococcus aureus.

BHI Broth 3ml in Screw Cap Tube in bulk

474CE.A



300 pieces 6 boxes of 50 pieces

eSwab®, liquid sample, swab sample

TSB™ (**Tryptic Soy Broth**) **Kit:** TSB™ is a generic multi-use enrichment broth without salts for aerobes, molds and yeasts.

TSB kit 2ml in Screw Cap Tube + 2 regular FLOQSwabs®



eSwab®, Fecalswab®, nose, throat or perineum sample

TSB broth 2ml in Screw Cap Tube in bulk

4U103N



300 pieces 6 boxes of 50 pieces eSwab®, Fecalswab®, liquid sample, swab sample

CAT Broth™: CAT broth™ is an enrichment broth for the isolation of Thricomonas vaginalis, Candida albicans, and other yeasts

CAT Broth 4ml Medium in 16x100 mm Screw Cap Tube + fiber swab in pouch

092CN



300 pieces 6 boxes of 50 pieces

swab sample



Thiol Broth™: Thiol medium™ is a universal enrichment broth for anaerobic bacteria and obligatory microorganisms. Its tube is also designed to enable turbidity check for strict anaerobes after 24h.



TSB Salt Enrichment Broth™: TSB Salt™ is a medium specifically intended for the isolation of S. aureus spp. After incubation, TSB Salt Broth™ is intended to enrich a specimen suitable for subculture on selective MRSA agar plates.



eMRSA Broth™: eMRSA broth™ is an enrichment broth for S. aureus spp, including MRSA and MSSA.

eMRSA™ Broth 2ml in Screw Cap Tube +1 regular FLOQSwabs® in pouch

086CE01.A



300 pieces 6 boxes of 50 pieces nose, throat or perineum sample

Selenite Broth™: Selenite Broth™ is a selective and enrichment broth for Salmonella enterica spp and Shigella sonnei. It is compatible both with stool or FecalSwab™ samples.

475CE.A

SELENITE Broth 2ml in Screw Cap Tube in bulk



300 pieces 6 boxes of 50 pieces

eSwab®, Fecalswab®

4E021N.A

SELENITE Broth 2ml in Rounded Bottom Tube



300 pieces 6 boxes of 50 pieces

eSwab®, Fecalswab®

GN Broth™: GN Broth™ is a selective and Enrichment broth for gram Negative bacteria, in particular Salmonella enterica and Shigella sonnei. It is compatible with untreated, eSwab® or FecalSwab™ samples.

085CU.A

 $\ensuremath{\mathbf{GN}}$ $\ensuremath{\mathbf{Broth}}$ 4ml in Screw Cap Tube in bulk



300 pieces 6 boxes of 50 pieces

eSwab®, Fecalswab®

LIM Enrichment Broth™: The LIM medium™ is a modification of Todd Hewitt Broth and is specifically intended for the isolation and enrichment of Group B Streptococcus.

LIM Broth 2ml in Screw Cap Tube

476CE.A



300 pieces 6 boxes of 50 pieces

 $eSwab^{\tiny{\circledR}}$



All above products are suitable for WASP® automation system. Please contact info@copangroup.com for additional information. *Suggested table. Please refer to your GLP procedures to choose the most appropriate device for the specific sampling site





This document may contain product information otherwise not accessible or valid in your country. Please be aware that Copan Italia S.p.A. does take any responsibility for accessing such information which may not comply with any valid legal process, regulation, registration or usage in the country of your origin. Product clearance and availability restrictions may apply in some Countries. Please refer to Copan website (www.copangroup.com) to view and/or download the most recent version of the brochure. This document is mainly intended for marketing purposes, always consult product insert for complete information. The use of this product in association with diagnostic kits or instrumentation should be internally validated by the user. ©2021 Copan Italia. All rights reserved. The trademarks mentioned herein are property of Copan Italia S.p.A.

Code: JMKPF017Roo.EN

