

Obturator

Product Information Card



for Surflo®/Versatus™ -W, Surshield® Versatus™ -S/W
and Surflo®/Versatus™ S-Plus

Description:

A sterile single use, closure device to be used in combination with the Surflo® -W, Versatus™-W, Surshield® Versatus™-S, Surshield® Versatus™-W, Surflo® S-Plus and Versatus™ S-Plus, I.V. catheter without injection port.

For plugging the peripheral catheter to prevent clot formation in the catheter during intermittent infusion or aspiration.

An obturator with the same colour as the catheter is available for each size except the 26G.

Does not contain components made of PVC or natural rubber latex.

Components:

Description	Raw Materials
Obturator	Polypropylene
Protector	Polyethylene

Product Range:



Product Code	Colour	Obturator Size		Length
SO+DM2419W	Yellow	24G	0.74 mm	19 mm
SO+DM2225W	Blue	22G	0.90 mm	25 mm
SO+DM2032W	Pink	20G	1.00 mm	32 mm
SO+DM1832W	Green	18G	1.30 mm	32 mm
SO+DM1845W	Green	18G	1.30 mm	45 mm
SO+DM1745W	White	17G	1.50 mm	45 mm
SO+DM1645W	Grey	16G	1.75 mm	45 mm
SO+DM1445W	Orange	14G	2.00 mm	45 mm

Packaging:

Single unit: thermoformed blister pack
Unit box: 50 pieces
Shipping carton: 500 pieces
 Gross weight = 4.6 kg
 Dimensions = 720 x 200 x 190 mm

Sterilization:

With ethylene oxide to a Sterility Assurance Level of at least 10⁻⁶

Obturator

for Surflo[®]/Versatus[™]-W, Surshield[®] Versatus[™]-S/W, Surflo[®]/Versatus[™] S-Plus

Storage:

Shelf life of 5 years

* The batch number and expiry date are printed on the paper of the thermoformed blister pack

* Do not store at extreme temperature and humidity

Precautions for Use:

Check integrity of the unit package before use

Do not reuse - do not resterilize

Dispose of safely after single use to avoid risk of infection

Status of Regulation:

Class IIa, sterile product following MDD 93/42/EEC

Notified body:

TÜV Rheinland 0197

Manufacturer:

TERUMO EUROPE NV

Interleuvenlaan 40

B - 3001 Leuven

Belgium



Terumo Corporation
+81 3 3374 8111

Terumo Europe NV
+32 16 38 12 11

**Terumo
Medical Products EMEA**
+32 16 38 12 11

EMEA Sales Offices

**Terumo Europe NV
African Business Division**
+32 16 38 13 08

**Terumo Europe NV
Benelux Sales Division
Belgium**

+32 16 339 25 80
The Netherlands
0800 0220396

**Terumo Europe NV
Emerging Markets Division**
+32 16 38 12 11

Terumo Deutschland GmbH
+49 6196 80 230

Terumo Deutschland GmbH
Switzerland
+41 56 419 10 10

Terumo Europe España SL
+34 902 10 12 98

Terumo France S.A.S.
+33 130 96 13 00

Terumo Italia S.R.L.
+39 0651 96 14 20

Terumo Russia LLC
+7 495 988 4740

Terumo Sweden AB
+46 3174 85 880

**Terumo Sweden AB
Denmark**
+45 7020 93 80

Terumo Middle East FZE
+971 4 292 0200

Terumo UK Ltd
+44 1276 480 440

**Terumo BCT Tibbi Cihazlar Dagim
ve Hizmetleri A.S.**
+90 216 590 0074

