

Terumo® 3-part Syringe



Product Information Card

Description:

Sterile 3-part syringe for single use of *Terumo Corporation Japan origin*.

The syringe is intended for the injection of fluids - immediately after filling - into the human body, for aspiration of fluids or for the sampling of blood or other body fluids.

Components:

Description	Raw Materials
Barrel	Polypropylene
Plunger	Polypropylene
Gasket	Thermoplastic elastomer
Lubricant	Silicone oil



Product Range:

Product Code	Size	Tip
SS*02LE1	2.5 ml	Concentric, Luer Lock
SS*02SE1	2.5 ml	Concentric, Luer Slip
SS*05LE1	5 ml	Concentric, Luer Lock
SS*05SE1	5 ml	Concentric, Luer Slip
SS*10LE1	10 ml	Concentric, Luer Lock
SS*30LE1	30 ml	Concentric, Luer Lock
SS*30ESE1	30 ml	Eccentric, Luer Slip

Technical Data:

Product Code	Nominal Capacity	Graduated Capacity	Increments
SS*02LE1	2.5 ml	2.5 ml	0.1
SS*02SE1	2.5 ml	2.5 ml	0.1
SS*05LE1	5 ml	6 ml	0.2
SS*05SE1	5 ml	6 ml	0.2
SS*10LE1	10 ml	12 ml	0.2
SS*30LE1	30 ml	35 ml	1.0
SS*30ESE1	30 ml	35 ml	1.0

Terumo® 3-Part Syringe

Packaging:

Single unit: blister pack (film/film)

Product Code	Unit Box	Shipping Carton		
		Quantity	Gross Weight	Dimensions (LxWxH)
SS*02LE1	100	2400	6.00 kg	578 x 386 x 436 mm
SS*02SE1	100	2400	9.44 kg	578 x 386 x 436 mm
SS*05LE1	100	1600	6.25 kg	544 x 404 x 444 mm
SS*05SE1	100	1600	10.30 kg	404 x 574 x 444 mm
SS*10LE1	100	800	15.00 kg	580 x 350 x 377 mm
SS*30LE1	50	400	16.00 kg	580 x 398 x 418 mm
SS*30ESE1	50	400	16.00 kg	580 x 398 x 418 mm

Sterilization:

With e-beam radiation to a Sterility Assurance Level of at least 10⁻⁶

Storage:

Shelf life of 5 years

* The batch number and expiry date are embossed in the film of the blister pack

* Do not store at extreme temperature and humidity

Precautions for Use:

Check integrity of the unit package before use

For single use only

Dispose of safely after single use to avoid risk of infection

Status of Regulation:

Class IIa, sterile product following MDD 93/42/EEC

Notified body:

TÜV Rheinland 0197

Manufacturer:

TERUMO CORPORATION JAPAN

EC representative:

TERUMO EUROPE N.V.

Interleuvenlaan 40

B - 3001 Leuven

Belgium



Terumo Corporation

44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

Laboratoires Terumo France S.A.

Parc Ariane - Bât. Uranus
1, rue Hélène Boucher
78284 Guyancourt Cedex, France
Tel.: +33 130 96 13 00
Fax: +33 130 43 60 85

Terumo Italia SRL

Via Simone Martini, 143/145
00142 Roma RM,
Italy
Tel.: +39 0651 96 14 20
Fax: +39 0650 30 407

Terumo Europe N.V.

Moscow Representative Office
BC Northern Tower
13Fr., 10 Testovskaya Street
Moscow 123317, Russia
Tel.: + 7 495 988 4740
Fax: + 7 495 988 4739

European Headquarters

Terumo Europe N.V.

Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

European Offices

Terumo Europe N.V.

Benelux Sales Division
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 39 25 80
Fax: +32 16 39 25 99
The Netherlands:
Tel.: 0800 0220396
Fax: 0800 0220414

Terumo Deutschland GmbH

Ludwig-Erhard-Strasse 6
65760 Eschborn, Germany
Tel.: +49 6196 80 230
Fax: +49 6196 80 23 200

Terumo Deutschland GmbH

Zweigniederlassung Spreitenbach
Bodenackerstrasse 3
8957 Spreitenbach, Switzerland
Tel.: +41 56 419 10 10
Fax: +41 56 419 10 11

Terumo Sweden AB

Nya Varvet, Byggnad 90
Sven Källfelts Gata 18
426 71 Västra Frölunda,
Sweden
Tel.: +46 3174 85 880
Fax: +46 3174 85 890

Terumo Denmark

Filial of Terumo Sweden AB
Tel.: +45 7020 93 80
Fax: +45 7020 94 80

Terumo Europe España SL

Avda. Juan Carlos I, Nº 13 - 7ª Planta
Edificio Torre La Garena
28806 Alcalá de Henares (Madrid),
Spain
Tel.: +34 9021 01 298
Fax: +34 9021 01 358

Terumo UK Ltd

Tamesis
The Causeway
Egham, Surrey
TW20 9AW, United Kingdom
Tel.: +44 1784 476 220
Fax: +44 1784 476 201

Terumo Corporation Dubai Branch

(Middle East and Africa)
Al Masraf Tower, 22nd Floor
P. O. Box 20291 Dubai
U.A.E.
Tel.: + 971 4 221 2220
Fax: + 971 4 221 3330

**Addresses subject to change,
please consult our Website.**

© Registered Trademark
Published by Terumo Europe N.V.
HO 19GB-0307TE-V(06.12)E