

RollBag[®]

Sterile sampling bag

Roll

Stick



and
that's it!


interscience

RollBag® is a sampling bag with a sticker closure for the collection of any liquid or solid sample, for the food industry, medical, animal health or environmental industries.

Developed to guarantee optimal sampling conditions, it preserves the integrity of the samples. RollBag® can be treated as soft waste which can be crushed and incinerated which allows the waste category to be downgraded to a more economical category.

Waterproof closure, with just a sticker!

Simply close your bag with the built-in sticker for a waterproof extra-resistant closure, without any extra tools.



1 Tear the pre-cut opening

2 Place the sample in the bag

3 Open the sticker & roll the bag 3 times

4 Close the bag hermetically

Lower your costs

RollBag® is metal free, meaning that it is **easier to eliminate than metallic wire bags.**

Indeed, the metallic wire can puncture the bags during autoclaving cycles and implies a disposal in the infectious risk biohazard waste stream in the category of sharp objects waste. This implies a high waste treatment cost.

Metal free RollBag® means you can autoclave it easily, or dispose of in a soft non-perforating waste biohazard class, that can be autoclaved, or incinerated and grinded. This allows the waste category to be downgraded to a more economical category!




Space saving

Rather than storing a large number of sampling containers, the use of RollBag® reduces storage and shipping costs.



500
RollBag®



500
plastic container



Water-tight closure
Easy & strong closing with the sticker



Resistant
Made in PolySilk®
Puncture proof



Extra strong welding
Resistant to blending



Safe
Gamma ray sterilization
Food contact approved



Pre-cut opening
Single sterile packaging



Marking zone
Easy identification of the sample



11 models for all applications

RollBag®													
Model	7.5 x 20	11 x 20	11 x 25	14 x 23	11 x 30	14 x 33	17.5 x 30	17.5 x 40	25 x 38	25 x 50	38 x 56		
Reference	145 072	145 112	145 115	145 142	145 113	145 143	145 173	145 174	145 253	145 255	145 385		
Ideal	For various types of any solid or liquid samples from the sampling site to the laboratory												
Width	cm	7.5	11	11	14	11	14	17.5	17.5	25	25	38	
	inch	3	4.3	4.3	5.5	4.3	5.5	7	7	10	10	15	
Height	cm	20	20	25	23	30	33	30	40	38	50	56	
	inch	8	8	10	9	12	13	12	16	15	20	22	
Thickness	µm	70	70	70	70	70	70	70	70	110	110	110	
	mil	3	3	3	3	3	3	3	3	4	4	4	
Max volume	mL	110	200	400	400	620	850	1100	1800	2500	4500	10000	
	oz	4	7	14	14	21	29	38	61	85	152	340	
Max blending volume	mL	65	160	180	200	145	350	650	550	1000	2500	5000	
	oz	2	5	6	7	5	12	22	19	34	84	170	
Optimal blending volume	mL	5-40	5-80	5-110	5-130	5-100	50-250	50-400	50-350	200-600	400-1000	500-3000	
	oz	0.15-1.5	0.15-3	0.15-4	0.15-5	0.15-4	1.5-9	1.5-14	1.5-12	7-21	13-34	17-102	
Compatibility with interscience lab blender	MiniMix®					MiniMix® & BagMixer®		BagMixer®		BagMixer® & JumboMix®		JumboMix®	
Nb of stickers	1	1	1	1	1	1	1	1	2	2	2		
Type of bag	Sampling bag												
Marking zone	✓												
Composition	PolySilk®: high grade HDPE (bag) - polypropylene (closure sticker)												
Gamma ray treated	5 to 12 kGy, with certificate												
Average shelf-life	5 years												
Temperature resistance	From -20°C to 70°C / -4°F to 158°F												
In compliance with	ISO 7218, ISO 6887 and FDA BAM (Bacteriological Analytical Manual)												
Approved for food contact	Regulation (EC) N° 1935/2004												
Storage	In a temperate and dry place, away from humidity and light												
Box dimensions	15.5 x 25.5 x 13 cm				35 x 21 x 10.5 cm				44 x 27.5 x 10 cm				41 x 33.5 x 11.3 cm
Gross weight	1.3 kg	1.7 kg	2.2 kg	2.5 kg	2.4 kg	3.5 kg	3.7 kg	5 kg	4.9 kg	6.3 kg	10.4 kg		
Box of	500	500	500	500	500	500	500	500	250	250	250		

Ask for free samples!

Certified production

Products manufactured for INTERSCIENCE by Interlab, an ISO 9001 certified company.



Your local distributor

interscience

PARIS

Phone: +33 (0)1 34 62 62 61 - Email: info@interscience.com

FRANKFURT

Phone: +49 611 7238 7770 - Email: sales.germany@interscience.com

BOSTON

Phone: +1 781 937 0007 - Email: sales.usa@intersciencelab.com

SHANGHAI

Phone: +86 (0)21-64739390 - Email: sales.china@interscience.com

SINGAPORE

Phone: +65 6977 7232 - Email: sales.asia@interscience.com

TOKYO

Phone: +81 3 6712 9715 - Email: sales.japan@interscience.com

www.interscience.com