



# TOPLACTINE LACTOSE BEVERAGE WITHOUT COLORANT

High concentrated lactose  
solution



A ready to use lactose solution,  
without colorant



The TOPLACTINE solution is used in lactose  
tolerance testing to diagnose **lactose intolerance**



◆ A solution for **Oral Lactose Tolerance Test**

- ready to use, easy to take
- standardized and accurate lactose concentration: each bottle of 300ml contains 50g of lactose
- pleasant orange flavour, **no colorant added**

There are two types of lactose tolerance tests: these tests must be performed under medical supervision

- ◆ **Hydrogen breath test.** This test measures the amount of hydrogen gas in your breath before and after you drink a **liquid containing lactose**. It is the most common way to test for lactose intolerance.
- ◆ **Glucose blood test.** This test involves a series of blood tests that measure the level of glucose (a type of sugar) in your blood before and after you drink a **liquid containing lactose**.



*Lactose intolerance is the inability to digest significant amounts of lactose, the primary sugar in milk. This inability results from a shortage of the enzyme lactase, which is normally produced by the cells that line the small intestine. Lactase breaks down lactose into simpler forms that can then be absorbed into the bloodstream during the digestive process. Common symptoms of lactose intolerance include nausea, cramps, bloating, gas, and diarrhea.*

Reference	Flavour	Lactose dose (g)	Bottle volume (ml)	Minimum Order Quantity (MOQ)
TOPLACTINE	Orange	50	300	20