

TOPSTAR GLUCOSE DRINK WITHOUT COLORANT

High concentrated glucose solution



A ready to use glucose solution, gluten-free and without colorant

TOPSTAR solutions allow diagnosis of:

- · Gestational diabetes
- Type 1 diabetes or diabetes mellitus,
- Type 2 diabetes
- Solutions for Oral Glucose Tolerance Test (OGTT)
 - o ready to use, easy to take
 - standardized and accurate glucose concentration
 - o available in 3 flavours: raspeberry, orange and lemon, without colorant



- OGTT measures the body's ability to metabolize glucose. It should be carried out by specialized personal, early in the morning and after 12 hours fasting. In these conditions, patient should drink a solution containing a known amount of glucose, in a period no longer than 10 minutes
- The patient's blood is tested in fasting and again 30 min, 1 hour, 2 hours and 3 hours after taking the high glucose drink.







Reference	Flavour	Glucose dose (g)	Bottle volume (ml)	Minimum Order Quantity (MOQ)
TOPSTARFRA75SC	raspeberry	75	200	30
TOPSTARCIT75SC	lemon	75	200	30
TOPSTAROR75SC	orange	75	200	30

